




POM Wireless Converter
Quickstart Guide

GLN2041

Setting



- The POM Wireless Converter works with Gleen sensors must be paired with POM Wireless Integrated Alarm to function correctly.
- 2. Make sure the sensor's cable connector is correctly inserted into the POM Wireless Converter socket.






Seat Belt or Sensor Pad

*Available Separately

Pairing

2



Pairing
Blue Button

| Continue |

- 1. Turn the GLN0111 POM Wireless Integrated Alarm and GLN2041 POM Wireless Converter ON.
- 2. Make sure Wireless Integrated Alarm and Wireless Converter are rapidly flashing green at the same time, before you press the Pairing buttons. Press the Blue button on the back of Wireless Integrated Alarm and press the Black button on Wireless Converter until the montior indicator light rapidly flashes green as well.
- When pairing is successful, the light on both devices should slowly flash green 3 times.
- 4. If pairing is unsuccessful, red light will flash slowly 3 times, turn all devices off and repeat from step 1 to step 3.
- Test the System by applying pressure on the sensor pad for 3 seconds, alarm will beep for 3 times meaning it is activated. Release pressure, the alarm will sound now.

FCC Notice

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

NOTE 1: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- -Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- -Consult the dealer or an experienced radio/TV technician for help.

NOTE 2: Any changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.