

The product description

Foreword

Thank you for using HEJIAJANKANG-SmartPillbox. Please read this manual carefully to fully understand how to use this device. I wish you happy! The company reserves the right to make changes to the contents of this manual without any prior notice.

Product parameters

Name: Box Home Health Network:

GSM wireless data terminal model: YH2

Weight: 129g

Battery: Lithium polymer battery 3.8V 1000mAh

Charging input: DC 5V 500mA

Working environment temperature: 0-45°C

Relative humidity: non-condensing state $60 \pm 25\%$ R.H

Precautions

1. Please take the medicine in the medicine cabinet with flashing light.
2. Keep the warehouse clean and dry, and clean it when changing medicines
3. After taking the medicine, please insert the medicine tank into the corresponding numbered grid and store the medicine in the dark
4. The kit is a medication assistive device. During the setting process, you should pay attention to the correct setting time of the alarm clock to ensure accurate medication.

Special Note: This product is not waterproof, please put it in a cool dry place, and avoid direct sunlight.

Before using this product, please read this chapter carefully, follow these instructions, and pay attention to these warnings. (1) The use of this product in the following places and weather is strictly prohibited:

1. Flammable and explosive places. Gas stations, fuel depots, chemical plants, places where blasting operations are being performed; areas where air contains particles, powders, metal powders and other chemicals or particles are potentially hazardous to explosions (such as the transportation or storage of tanks, fuel, or chemicals) facility).

2. Aircraft cockpit. Please turn off this product when taking off the aircraft. To prevent interference with the aircraft communication system, it is forbidden to use this product during flight.

3. Medical products and public medical facilities. Any radio transmitting device (including this product) may interfere with the function of improperly protected medical devices. If there is a problem or if it is uncertain whether the medical device can shield external RF energy, consult your doctor or medical device manufacturer.

4. Other places that prohibit the use of unlimited transmission capabilities of electronic devices.

5. Thunderstorm weather. Thunderstorm weather may cause malfunction or shock hazard.

(b) Keep away from the following interfering medical equipment:

1. Pacemakers. When opening the product, keep the distance between this product and the pacemaker greater than 15.3 cm (approximately 6 inches) at all times. Do not place this product in a jacket pocket. Use this product on the other side of the pacemaker to reduce potential interference. If you suspect that interference has occurred, turn off the product immediately.

2. Hearing aids and implants into the cochlea. Some digital wireless devices may cause some interference with certain hearing aids or implanted cochleas. If such interference occurs, consult an authorized person and refer to relevant medical standards.

3. Other medical equipment. If you and other people around you also use other medical devices, consult an authorized medical staff or equipment manufacturer to determine if they can adequately shield the surrounding RF signals.

(III) Other safety precautions

1. Do not attempt to open the main unit of this product. If handled or handled improperly, this product and its contents may affect human health or pollute the environment.

2. Do not use this product in dusty, wet, dirty, or near magnetic fields, as this may cause internal circuit failure of the product.

3. The use of radio frequency equipment including this product may cause some interference to electronic equipment with poor shielding performance. Consult the manufacturer of the equipment before using the product if necessary.

4. Keep all parts such as chargers and data cables out of reach of children and pets to prevent suffocation or other dangerous situations when children and pets swallow parts.

5. The normal operating temperature range of this device is -10 degrees Celsius to 45 degrees Celsius. Outside this temperature range may cause product malfunction, fire or explosion.

6. Do not allow children to charge the cartridge or touch the charging socket on their own to prevent the risk of electric shock.

7. Do not drop or bump the charger.

The company and its manufacturers are not responsible for any consequences of improper use or non-compliance with the product used in this manual. The company reserves the right of final interpretation of this specification to the extent permitted by law.

FCC Warning

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions

(1) this device may not cause harmful interference, and
(2) this device must accept any interference received, including interference that may cause undesired operation. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. Note 1: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC RF Exposure Information and Statement:

H947 has been tested against this SAR limit. This device was tested for typical body-worn operations with the back of the product kept 0mm from the body. To maintain compliance with RF exposure requirements, use accessories that maintain a 0mm separation distance at the user's body . The use of belt clips, holsters and similar accessories should not contain metallic components in its assembly. The use of accessories that do not satisfy these requirements may not comply with RF exposure requirements, and should be avoided. Body-worn

Operation This device was tested for typical body-worn operations. To comply with RF exposure requirements, a minimum separation distance of 0mm must be maintained at the user' s body including the antenna. Third-party belt-clips, holsters, and similar accessories used by this device should not contain any metallic components. Body-worn accessories that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna.