## English User Guide

Thank you for purchasing the CYCO product. This user guide outlines the operating procedures. Please read it carefully to make sure you are using your product properly and retain it for future reference.

### How to use

- A. CYCO Smart Pill Box
  1. Connect ❶ to ❷
  2. Connect ❸ to a USB port with a power supply (e.g. AC adaptor with USB power supply).
  3. The smart pill box's LED indicator will flash red while charging.
  4. When fully charged, the LED indicator will become flash green.
- B. Connecting the smart pill box via APP and Bluetooth
  1. Enable APP and Bluetooth on the device you want to connect the smart pill box with.
  2. Hold down ❹ for 3 secs until the LED indicator flash green. This will start the connecting device and the LED indicator will flash blue quickly.
  3. Once your device and the smart pill box are paired, your device will display "SmartPillbox", in the list of Bluetooth devices to connect to and the LED indicator will become flash blue.
  4. Under the installation guidelines from APP, you can build a healthy lifecycle easily.
- C. Other cases
  1. When the smart pill box's battery is low, the LED indicator will flash red.
  2. When the smart pill box is not being connected to APP and Bluetooth, the LED indicator will flash green.
  3. When the smart pill box is connected to APP and Bluetooth, the LED indicator will flash blue.
  4. Press ❺ for mute function when "beep" alert.
  5. Press ❻ for 10 secs until the smart pill box is off.

批注 [L1]: With APP

Note: Please disconnect any previous device that was paired with the smart pill box. If a previous device is still paired via Bluetooth and is within range, you will not be able to connect another device until it has been unpaired.

### Additional Information

Please register your CYCO product online at [www.rightcyco.com](http://www.rightcyco.com).

Please visit our website [www.rightcyco.com](http://www.rightcyco.com) or contact our customer service center at [cs@qualife.com.hk](mailto:cs@qualife.com.hk) to learn more about our new products and special offers.

#### FCC Statement

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment

#### IC WARNING

This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

- (1) This device may not cause interference.
- (2) This device must accept any interference, including interference that may cause undesired operation of the device.

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