

WARNINGS

BEFORE USING OR CHARGING THE ELITE, READ ALL THE INSTRUCTIONS AND CAUTIONARY MARKINGS IN THIS MANUAL, ON THE CHARGER, AND ON THE ELITE. The ELITE is intended for commercial or home use. We strive to make the ELITE as safe as possible. This is an advanced mechanical device with electric components. If the ELITE and its accessories are not used or maintained properly, there is a risk of fire, electric shock, or injury. Do not place your finger or any objects near the metal piston above the attachment while the ELITE is in use or immediately after use. The metal piston can become extremely hot. Do not use any attachments other than the ones provided by Theragun and only as instructed by Theragun. Do not use the ELITE on your head or near your genitals. Do not use the ELITE above the waist during pregnancy. Do not get hair caught in any moving part of the ELITE. Do not use the ELITE above your Adam's apple or C4. Do not get the ELITE wet. Only clean the ELITE with a magic eraser or a lightly damp towel. Do not block the vents of the motor. Never attempt to disassemble the ELITE. Never store the ELITE at temperatures above 40°C/104°F. Only charge the ELITE when the ambient temperature is between 0°C/32°F and 40°C/104°F. Never expose the ELITE or its charger to water. For long-term storage, the ELITE should be fully charged. This product is not intended to diagnose, treat, or prevent any diseases. When using the ELITE, these basic precautions should always be followed including:

1. USE ONLY AS INSTRUCTED. Use the ELITE as described in this ELITE Instruction Manual. Use only Theragun recommended accessories and replacement parts. Do not carry out any maintenance other than what is shown in the ELITE Instruction Manual or advised by Theragun.
2. NOT FOR CHILDREN. The ELITE and charger are not intended for use by young children or persons with reduced physical, sensory or reasoning capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction by a responsible person. Do not allow the ELITE to be used as a toy. Close attention is necessary when used by or near children. Children should be supervised to ensure that they do not play with the ELITE or charger.
3. REDUCING RISK OF INJURY. Charge the ELITE only with the Theragun ELITE charger. Other chargers may cause personal injury or damage. Do not connect the ELITE to a different power supply plug or car cigarette lighter. The internal battery will be permanently disabled or damaged. Only use the appropriate Theragun Wireless Charger to wirelessly charge the ELITE.
4. CHARGING LOCATIONS. The ELITE should be charged indoors in a well ventilated, dry location. Do not charge the ELITE outdoors, in a bathroom or within 10 feet (3.1 metres) of a bathtub, or pool. Do not use the ELITE or charger on wet surfaces and do not expose the charger to moisture, rain, or snow. Do not use the ELITE or charger in the presence of explosive atmospheres (gaseous fumes, dust, or flammable materials). Sparks may be generated, possibly causing a fire.
5. UNPLUGGING THE CHARGER. Pull the plug, not the cable, to reduce the risk of damage to the electrical plug and cable. Never carry the charger by its cable. Keep the cable from heat, oil, and sharp edges. Do not use or charge with a damaged cable. Do not stretch the charger cable or place the cable under strain. Keep cable away from heated surfaces. Do not handle the charger, including the charger terminals or ELITE with wet hands.
6. CHARGER CARE. Unplug the charger when not in use. Unplug the charger from the outlet when not in use for extended periods. Make sure the charger cable is located so that it will not be stepped on, tripped over, or otherwise subjected to damage or stress. Have a damaged charger replaced immediately.
7. DO NOT OVERCHARGE. Do not leave the ELITE connected to the charger for more than an hour after the ELITE has been fully charged. The ELITE includes a protection system to avoid the risk of overcharging. However, overcharging the ELITE may reduce its life over time. Only charge with AC Adaptor: EA1046SAR or with the Theragun ELITE Wireless Charger Stand.
8. DO NOT BLOCK MOTOR VENTS. Do not put any obstructions into vents. Do not use if any vent is blocked. Keep free of dust, lint, hair, and anything that may reduce airflow. Keep hair, loose clothing, and fingers from vents and moving parts.
9. DO NOT BURN OR INCINERATE THE ELITE. The internal battery may explode, causing personal injury or damage. Toxic fumes and materials are created when battery is burned.
10. DO NOT CRUSH, DROP, OR DAMAGE THE ELITE CHARGER. Do not use a charger that has received a sharp blow, been dropped, run over, or damaged in any way.
11. BATTERY CHEMICALS CAUSE SERIOUS BURNS. Never allow the internal battery to come into contact with the skin, eyes, or mouth. If a damaged battery leaks chemicals, use rubber or neoprene gloves to dispose of it. If skin is exposed to battery fluids, wash with soap and water, then rinse with vinegar. If eyes are exposed to battery chemicals, immediately rinse with water for 20 minutes and seek medical attention. Remove and dispose of contaminated clothing.
12. DO NOT OPERATE UNDER BLANKET AND PILLOW. Excessive heating can occur and cause fire, electric shock, or injury to persons.
13. STORING THE ELITE AND CHARGER. Store in a cool, dry place. Do not store the ELITE where temperatures may exceed 40°C/104°F such as in direct sunlight, in a vehicle, or metal building during the summer.
14. DISPOSING OF BATTERIES. The internal Theragun Lithium-ion batteries for the ELITE are more environmentally friendly than some other types of batteries (e.g., nickel-cadmium). Always dispose of the ELITE according to government, regional, and local regulations. Contact a recycling agency in your area for

recycling locations. Even discharged batteries contain some energy. Before disposing, use electrical tape to cover the ELITE's 12V connector to prevent the battery from shorting, which could cause a fire or explosion.

15. DO NOT DISASSEMBLE. Disassembly or incorrect reassembly may result in the risk of electric shock, fire, or exposure to battery chemicals. The warranty will be void if the ELITE or charger are disassembled or if any parts have been removed.
16. SERVICE. If the ELITE or charger is not working properly, has received a sharp blow, has been dropped, damaged, left outdoors, or dropped into water, do not use and contact Theragun by emailing info@theragun.com for repair or service. Do not attempt to repair or disassemble the Theragun ELITE which may result in an electric shock or fire.
17. DO NOT USE WHILE BATHING OR IN THE SHOWER, TUB, OR SINK. Do not place or store the ELITE where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid. Do not reach for an appliance that has fallen into water. Unplug immediately.
18. THERMAL LIMITER. The ELITE employs an automatic reset thermal limiter that shuts off the ELITE to prevent overheating and fire. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

LIMITED WARRANTY

The ELITE and charger have a limited 1-year warranty. To have the ELITE, a part repaired, or replaced under warranty, it must be returned with postage prepaid and insured, to Theragun. This warranty does not apply to damage that Theragun determines to be from repairs made or attempted by anyone other than Theragun authorised personnel, misuse, alterations, abuse, normal wear and tear, lack of maintenance, or accidents. For the full details of the warranty, please go to www.theragun.com/warranty. For any warranty issues, please email info@theragun.com.

IC

1. This device complies with Industry Canada RSS standard(s). Operation is subject to the following two conditions:
(1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.
2. Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment complies with ISED RF radiation exposure limits set forth for an uncontrolled environment.

1. Le présent appareil est conforme aux CNR d'Industrie Canada applicables aux appareils radio exempts de licence. L'exploitation est autorisée aux deux conditions suivantes:
(1) l'appareil ne doit pas produire de brouillage, et (2) l'utilisateur de l'appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d'en compromettre le fonctionnement.
2. Les changements ou modifications non expressément approuvés par la partie responsable de la conformité peuvent annuler le droit de l'utilisateur à utiliser l'équipement.

Cet appareil est conforme aux limites d'exposition de rayonnement RF ISED établies pour un environnement non contrôlé.