

---

## WARNINGS

BEFORE USING OR CHARGING THE PRO, READ THIS MANUAL AND THE CAUTIONARY MARKINGS ON THE CHARGING STAND AND THE PRO. CHARGE THE PRO BEFORE FIRST USE. We strive to make the PRO and Charging Stand as safe as possible. This is an advanced mechanical tool with electric components. If the PRO and its accessories are not used or maintained properly, there is a risk of fire, electric shock, or injury. When using the Charging Stand, these basic precautions should always be followed:

1. USE ONLY AS INSTRUCTED. Use the Charging Stand as described in this Instruction Manual. Use only Theragun recommended accessories and replacement parts. Do not carry out any maintenance other than as advised by Theragun.
2. NOT FOR CHILDREN. The PRO and Charging Stand are not intended for use by young children or persons with reduced physical, sensory or reasoning capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction by a responsible person. Close attention is necessary when used by or near children. Children should be supervised to ensure that they do not play with the PRO or the Charging Stand.
3. REDUCING RISK OF INJURY. Charge the PRO only in an authorized Theragun charger. Other chargers may cause personal injury or damage. Do not connect the PRO to a power supply plug or car cigarette lighter. The battery will be permanently disabled or damaged.
4. CHARGING LOCATIONS. The PRO should be charged indoors in a well ventilated, dry location. Do not charge the PRO outdoors, in a bathroom or within 10 feet (3.1 meters) of a bathtub, or pool. Do not use the PRO or Charging Stand on wet surfaces and do not expose the charger to moisture, rain, or snow. Do not use the PRO or Charging Stand in the presence of explosive atmospheres (gaseous fumes, dust, or flammable materials). Sparks may be generated, possibly causing a fire.
5. UNPLUGGING THE CHARGING STAND. Pull the plug, not the cord, to reduce the risk of damage to the electrical plug and cord. Never carry the Charging Stand by its cord. Keep the cord from heat, oil, and sharp edges. Do not use or charge with a damaged cord. Do not stretch the charger cord or place the cord under strain. Keep cord away from heated surfaces. Do not handle the Charging Stand, including the charger terminals with wet hands.
6. CHARGING STAND CARE. Unplug the Charging Stand when not in use. Unplug the Charging Stand from the outlet when not in use for extended periods. Make sure the charger cord is located so that it will not be stepped on, tripped over, or otherwise subjected to damage or stress. Have a damaged charger replaced immediately.
7. DO NOT BURN OR INCINERATE THE CHARGING STAND OR THE PRO.
8. DO NOT CRUSH, DROP, OR DAMAGE THE CHARGING STAND. Do not use a charger that has received a sharp blow, been dropped, run over, or damaged in any way.
9. STORING THE CHARGING STAND. Store in a cool, dry place. Do not store the Charging Stand where temperatures may exceed 40°C/104°F such as in direct sunlight, in a vehicle, or metal building during the summer.
10. DO NOT DISASSEMBLE. Disassembly or incorrect reassembly may result in the risk of electric shock, fire, or exposure to battery chemicals. The warranty will be void if the Charging Stand is disassembled or if any parts have been removed.
11. SERVICE. If the Charging Stand is not working properly, has received a sharp blow, has been dropped, damaged, left outdoors, or dropped into water, do not use and contact Theragun by emailing [info@theragun.com](mailto:info@theragun.com) for repair or service. Do not attempt to repair or disassemble it which may result in an electric shock or fire.
12. DO NOT USE WHILE BATHING OR IN THE SHOWER, TUB, OR SINK. Do not charge the PRO where it can fall or be pulled into a tub or sink. Do not place in or drop into water or other liquid. Do not reach for an appliance that has fallen into water. Unplug immediately.
13. The PRO can not be used while charging.

### LIMITED WARRANTY

The Charging Stand has a limited one-year warranty. For any warranty issues or questions, please email [info@theragun.com](mailto:info@theragun.com).

---

changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

This equipment should be installed and operated with minimum distance 20cm between the radiator and your body