

80mm

150mm

# WYZE SCALE


## QUICK START GUIDE

**WELCOME TO WYZE**

Hi! Thanks for getting your hands on a Wyze Scale, and if this is your first Wyze product, welcome to the Wyze family! You just made a commitment to a better, more connected lifestyle. This guide will help you get things set up.


**IN THE BOX**

- Wyze Scale × 1
- 1.5V AAA Batteries × 4
- Quick Start Guide × 1




**SETUP**


- Install 4 included AAA batteries. Please note the positive and negative signs to ensure the batteries are installed correctly.




- Download the Wyze App and login or create an account.






- Follow the instructions on the Wyze App to add the scale. To connect, please enable your phone's Bluetooth, tap the scale to activate it, and have your phone nearby.



- After the connection is established, place the scale on a hard, flat surface and step onto the scale with bare feet to begin the measurement.












**TIPS**

-  Place on a hard, flat surface
-  Don't wear socks
-  Don't keep your knees together

**NOTE**

- You can share the Wyze Scale with other family members by following the sharing instructions inside the Wyze App.
- Your scale may need to be calibrated to show accurate weight measurement. Place the scale on a hard, flat surface and re-calibrate it (step on and then step off, allowing the scale to return to "0.0 lb").
- If you work out more than 8 hours a week and have a resting heart rate below 60 bpm, you may need to use the Athlete mode.
- Make sure your age, sex, and height entered during the setup process are correct to get the most accurate results.
- Open the Wyze App and connect to the Wyze Scale before taking measurements to ensure the data is logged on the correct profile.
- If the Wyze Scale is shared by users with similar weights, make sure the scale is connected to the Wyze App of the user taking the measurement or data may not be recorded on the correct profile.
- People that are under 18 or over 65 years of age, that are pregnant, that have osteoporosis or very low bone density, that have edema, or that have artificial dialysis should use the Wyze Scale only as a reference as the body fat rate could be inaccurate due to deviations from standard averages.

**LED DISPLAY**

-  Suspected error. Please step on the scale bare feet and clean the scale with damp cloth if needed.
-  The batteries are low. Please replace the batteries.
-  The batteries are critically low. Please replace the batteries.
-  Module connection error. Please reinstall the batteries.
-  Weight not in supported range. Please measure again.
-  Unstable measurement. Please measure again.
-  Hardware error. Please reinstall the batteries and contact customer support if needed.
-  Firmware update error. Please try again.
-  App connection error. Please restart the App and try again.

**SPECIFICATIONS**

- Product Name: Wyze Scale
- Model: WHSCL1
- Product Size: 306x306x26.6mm
- Product Weight: 4.07lb/1.85kg
- Power: 4 X 1.5V AAA batteries (included)
- Display Units: lb/kg
- Compatibility: Android 5.0+, iOS 9.0+
- Measurement Range: Weight: 11 - 400 lb/5 - 180 kg
- Operating Temperature: 41°F - 104°F (5°C - 40°C)
- Operating Humidity: 20 - 85%RH
- Storage and Transport temperature: -4°F - 131°F (-20°C - 55°C)
- Storage and Transport humidity: ≤90%RH

**CAUTION**

Never use this product in combination with medical electronic devices such as:

- Medical electronic implants such as pacemakers.
- Electronic life support systems such as artificial heart/lungs.
- Portable electronic medical devices such as electrocardiographs.

Pregnant women need to consult their healthcare provider before use. Some physical conditions could affect hydration levels that may lead to inaccurate results. Please consult your healthcare provider for more information.

**GENERAL SAFETY & PRECAUTIONS**

- Do not stand on the edge of scale, as you may fall over or influence the accuracy of the measurement.
- Do not use the scale on a tiled or wet floor, as this may result in a fall.
- Make sure the surface of the scale is clean before you step onto it, to prevent any slippage.
- Treat your scale with care. Do not drop it or jump on it. The scale is designed to be stood on; misuse or abuse may break the electronic sensors, cause injury, or result in inaccurate measurements.
- Avoid using this scale near strong magnetic fields, such as microwave ovens, etc.
- This device complies with part 15 of the FCC Rules. Its operation is subject to the following two conditions:
  - this device may not cause harmful interference, and
  - this device must accept any interference received, including interference that may cause undesired operation.
- Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

—Reorient or relocate the receiving antenna.

—Increase the separation between the equipment and receiver.

—Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.

—Consult the dealer or an experienced radio/TV technician for help.

8. This device contains licence-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's licence-exempt RSS(s). Operation is subject to the following two conditions:

(1) This device may not cause interference. (2) This device must accept any interference, including interference that may cause undesired operation of the device.

L' émetteur/récepteur exempt de licence contenu dans le présent appareil est conforme aux CNR d' Innovation, Sciences et Développement économique Canada applicables aux appareils radio exempts de licence. L' exploitation est autorisée aux deux conditions suivantes :

(1) L' appareil ne doit pas produire de brouillage; (2) L' appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d' en compromettre le fonctionnement.

9. This equipment complies with IC RSS-102 radiation exposure limits set forth for an uncontrolled environment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

Cet équipement est conforme aux limites d'exposition aux radiations IC CNR - 102 établies pour un environnement non contrôlé. Cet émetteur ne doit pas être situé au même endroit ou fonctionner avec une autre antenne ou un autre émetteur.

10. Bluetooth® associated logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Wyze Lab Inc. is permitted under license.

Other trademarks and trade names are those of their respective owners.




**CARE & MAINTENANCE**

- Avoid high temperatures and irradiation.
- If the scale is stored close to freezing point, allow it to acclimate to room temperature before use.
- Remove the batteries if the scale is not going to be used for 3 months or more.
- Clean the scale with a soft damp cloth. Do not use abrasive or solvent-based cleaners, and never immerse it in water.
- The scale will maintain its safety and performance features for at least 10,000 measurements.

**WARRANTY INFORMATION**


The Wyze Scale is warranted to be free from defects in materials and workmanship appearing within 1 year from the date of purchase, when used in accordance with the instructions provided. The above warranties extend only to the original retail purchaser. We will, at our option, repair or replace without charge any product covered by the above warranties. Repair or replacement is our only responsibility and your only remedy under the above warranties.

**EXPLANATION OF SYMBOLS ON UNIT**

-  Warning.
-  Please recycle where facilities exist. Waste electrical products should not be disposed of with household waste. Check with your local authority or retailer for recycling advice.
-  Keep dry.

**NEED HELP ?**

Contact Our Support Team:  
[www.wyze.com/support](http://www.wyze.com/support)  
[support@wyze.com](mailto:support@wyze.com)  
 1-844-WYZECAM



Connect With Wyze  
<https://wyze.com/community>

图纸提供单位

型号

WHSCL1

名称

快速指南

版本号

V1.0

规格

WHSCL1-KSZP01 V1.0

设计时间

2019.09.26

设计师

李杰

材质工艺要求

材质要求: 105g无光铜版纸

制作工序: 印刷+过水性消光油 +模切  
 正背印刷  
 风琴折 折 5 下  
 成品尺寸: 80x150mm  
 展开尺寸: 480x150mm

工艺要求: 色相正确、印迹牢固、  
 套印准确、各色套印、  
 不露杂色, 套印误差 ≤ 0.1mm  
 模切走位 ≤ ± 0.5mm  
 本件要求通过RoHS测试

颜色及 专色: K

比例 1:1 单位 mm