




Instagram: [defi\\_smart\\_official](#)

E-mail: [gasian@gasian.co.kr](mailto:gasian@gasian.co.kr)



## DEFI Smart Pillbox L

Instruction

## What is included



Smart pillbox L



USB-C charging cable

## Information

Product name	Smart Pillbox L
Model	GPBL-2.0
Power	3.7V Lithium-ion battery
Bluetooth frequency	2.4GHz
Temperature range	-4°F (-20°C)

## Safety information & Key features


**Safety information** ① If the product is not turned on, charge it with the enclosed charger.  
② if the alarm does not go off when it's connecting the pillbox, please retry after turn off the application.

**Key features** ① DEFI Smart Pillbox helps people to take their medicine on time by their own alarm setting.  
② Set an alarm on DEFI application and it will connect with a smartphone.  
③ DEFI Smart Pillbox needs Bluetooth connection with your smartphone.  
④ It extends battery life up to 30 days on a single charge.

## Installing the DEFI app




① search and download "DEFI" on Google Play or APP store.



③ After connecting with your smart pillbox, set an alarm.



② Ensure that your mobile's Bluetooth is turned on for the app to connect with your DEFI Smart Pillbox and create your DEFI account.



④ Check the operation by filling up your pillbox (IR sensor) and history (open/closed).


## Warranty Form

Product name	Smart Pillbox L
Model	GPBL-2.0
Serial number	
Place of purchase	
Purchase date	

## Key features


### 1 Bluetooth

Connect to you pillbox via Bluetooth with our mobile app




### 2 Consumption and refill alarm

Pillbox and smartphone beep when it's time for your pills and send you a notification if the pillbox is empty.




### 5 Proximity and disconnected Alarm

The alarm goes off when the distance between the smartphone and smart pillbox is over 33ft (10m), letting you know that you may be forgetting to take the pillbox with you.




### 6 Overdose alert

If the pillbox is accessed too often within a short period according to the pillbox's usage history, you'll receive an alert to prevent accidental double dosing.




### 3 Consumption history

the pills' existence can be shown on the consumption history as taken, missed, or unknown. Also, the pillbox automatically records each time you access (open/close) the pillbox.



### 4 IR sensed reminder

The IR sensors in each divided space can send you a refill notification when the compartment is empty. (Note: IR sensing can't detect soft capsule pills).



## Federal Communication Commission Interference Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

- (1) This device may not cause harmful interference, and
- (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment. This transmitter must not be co-located or operating in conjunction with any other antenna or transmitter.

## Radiation Exposure Statement:

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with 0.5cm between the radiator & your body.