

## *Instructions for use*

### **Power on**

In the off state, hold down the power button for 1.5 second

### **Power Off**

Hold down the power button for 1.5 seconds

### **Control output**

Press the up key to increase the output voltage. Long press on the build, can quickly improve the output voltage,

Short press the key to reduce the output voltage. Long press build, can quickly reduce the output voltage,

LED prompt

When the machine does not output,LED light breathing display.

After the output of the machine is turned on, the LED light is always on

### **LED prompt**

When the machine does not output,LED light breathing display.

After the output of the machine is turned on, the LED light is always on.

## **Electric quantity**

The machine will shut down automatically when the battery is too low.

When the battery is too low, the system will prompt low power and then automatically shut down.

## ***Bluetooth matching instructions***

Step 1: put the machine and foot into the matching mode at the same time

When the machine is turned off, press three buttons on the machine for 1.5 seconds at the same time. The machine enters the matching mode. The Bluetooth icon flashes on the screen.

In the off state, press the pedal control button and the power button at the same time for 1.5 seconds, and the pedal enters the matching mode. The LED light alternates red and green and flashes rapidly.

Step 2: Wait for the match to complete

The matching of the machine to the foot is done automatically. When finished, the machine will display a success sign. The LED indicator on the foot no longer flashes rapidly.

Step 3: Boot automatic connection

When the match is successful

After the machine and foot start, it will be automatically connected.

When the connection is successful, the Bluetooth icon will appear on the machine screen. The LED indicator of the foot is on.

FCC Warning Statement: Changes or modifications to this unit not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC Radiation Exposure Statement

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.