


说明书尺寸: 长100宽80

Visual Water Flosser User's Guide



C2

PACKAGE LIST

- Body *1
- Visual nozzle *1
- Standard nozzle *2
- Large elbow nozzle *2
- Manual *1
- Charging Cable *1

Product Description

After extensive adaptability research and repeated testing, we have made more precise function settings for people with oral sensitivity. The original visualization design enables real-time observation of tooth conditions during tooth rinsing, bringing an intuitive use experience. The innovative knob gear design provides finer gears to adjust water pressure and water flow in FLOSS/MASSAGE mode. To meet the needs of more different users, users can control it as they please, further improving the user experience.

Taking into account both cleaning power and gum protection: Innovative elastic water flow tooth cleaning technology outputs high-frequency flowing cleaning power to deeply clean the blind areas of the mouth such as between teeth and around braces. It is also equipped with 2-second crescendo technology. After the power is turned on, the water flow will gradually increase from weak to strong. Reduce the impact on gums and better protect gums.

Redefining practicality and thoughtfulness: Subverting the capacity of the traditional vertical dental rinsing, it is newly upgraded to a 350ml super large water tank, which can clean the entire mouth with one injection of water. The newly upgraded indicator light design makes it easy to grasp the current mode and battery status.

Safety and Maintenance

Do not use without nozzle.
Please make sure you fully understand the usage and techniques of this product before use.
Do not press the nozzle hard against your gums or teeth to avoid damaging your gums.
This product is not suitable for children and infants. Please store this product out of the reach of children or infants and do not use it as a toy.

Instructions for use

Nozzle installation and removal:

•Install
Insert the nozzle into the hole in the fuselage. When you hear the popping sound of the plug-in button, it is installed in place.

•Disassembly
Press and hold the nozzle insertion and extraction button, gently pull out the nozzle, and store the nozzle.

▲ When the product is running, it is prohibited to press the nozzle insertion and extraction button.

Nozzle Function

- Visual nozzle**
Observe the cleaning status of the dental cavity in real time
- Standard nozzle**
Suitable for daily cleaning of interdental debris and soft dent
- Large elbow nozzle**
Suitable for cleaning residue and soft dirt on the inside of teeth

Install App

Step 1:
Download "SUNUO" App (Android 4.4, iOS 9.0 or above are supported)
For IOS: Download and install the App"SUNUO" from App Store.
For Android: Download and install the App"SUNUO" from Google Play Store.
Optional: Scan the QR Code below to download the App.

Step 2:
Short press * CONNECT SMART DEVICES *

Step 3:
Connect wifi name: "SUNUO-XXXX" successfully and get back to app and click the "Connected", enter into the viewfind interface

" " Detailed Explanation

Water injection method:

- Punch water filling cover
- Keep the product in a horizontal state and pour water into the water tank from the water injection port.
- After injecting an appropriate amount of water, close the water injection cap.

▲ It is recommended to use warm boiled water or tap water below 40°C. Never pour mouthwash directly into your mouth. Never inject chemical products such as toothpaste, salt, solid mouthwash tablets, etc., otherwise they may become clogged or damaged.

Features:

- Off button**
Press and hold for 2S to turn on/off
Short press to pause
- Mode key**
Switching between FLOSS/MASSAGE modes
- Power Indicator**
The green light is always on when the battery is sufficient, the red light is flashing when the battery is low, and the red light is always on when charging.
- FLOSS mode light**
When working in this mode, the green light is always on and flashes when paused.
- MASSAGE mode light**
When working in this mode, the green light is always on and flashes when paused.

6 Gear switch
Pull it clockwise to increase the gear, and pull it counterclockwise to decrease the gear.
There are 7 levels in FLOSS mode, with intensity options from 8G-9G-10G-11G-12G-13G-14G.
There are 5 gears in MASSAGES mode, with intensity options from 10G/6G-11G/6.6G-12G/7.2G-13G/7.8G-14G/8.4G
The gear switch is self-returning, with a 45-degree rebound swing in both clockwise and counterclockwise directions.
A beeping sound will be heard when increasing or decreasing gears. If it exceeds the minimum or maximum gear, there will be no "beep" sound when the switch is pulled out.

▲ **Kind tips**
1. Mode adjustment and gear adjustment, with memory function. The next time you use it, the machine will automatically jump to the gear it was in last time.
2. Due to the working characteristics of the pump, the product will stop working continuously for up to 2 minutes. You need to press the power button to start. If the pause time exceeds 15 minutes, the machine will automatically shut down.

Tips to use

- Before use, it is recommended to select the teeth cleaning mode first, then place the nozzle against the gap between the teeth, and then click the start button to start rinsing. It is recommended to keep the lips slightly closed when using it to avoid water splashing and damaging the gums.
- When using it for the first time, adjust the gear according to your personal needs. Pull it once to adjust the gear. The gear switch is a rebound rotation. You need to release your finger for the switch to reset.

Adjust gear

Cleaning and maintenance

Clean after use
After use, it is recommended to wash and dry it thoroughly before storing it to avoid the growth of bacteria.

- Stop the device from running
- Open clockwise

Clean after use

- Dentists recommend replacing nozzles every three months
- If the top of the nozzle is deformed, please replace it with a new one in time to avoid affecting use.

Straw:

- Do not bend or pull the spray and suction pipe.

Product body:

- When cleaning the main body of the product, make sure that the silicone plug of the charging port is installed in place throughout the process to prevent water from infiltrating.
- Stains on the surface of the product can be removed with a cloth dipped in neutral detergent, and then wiped clean with a dry cloth. If the product body is very dirty, it can be cleaned with water. Do not use diluent, gasoline, alcohol, etc. to wipe this product.
- Do not use hot water above 50°C for cleaning.

Basic parameters

Product number	Rated power	Input parameters
C2	5W	5V = 1A
Nozzle diameter	Pulse frequency	Water pressure
0.6mm	1500 times/minute	40~130PSI
Water tank capacity	Battery capacity	Noise
350ML	2000mAh	< 70dB
Charging port	Charging time	Waterproof level
Type-C	≈2H	IPX7
Dimensions (mm)	Operating mode	Net Weight
81.5*69*233mm	2MODEL/12gear	300G
Pixel	WiFi	Viewing angle
30W	2.4G	78°
Focusing range	Auxiliary lighting	Video resolution
20-100mm	6LEDs	1080P

Water intake:

- Please check the water inlet regularly to ensure that no stains remain, otherwise it may cause water leakage or breed bacteria and affect your health.

Water tank:

- Please inject clean water every time you use it. It is recommended that you unplug the water tank part and clean the inner wall of the water tank after each use. This will reduce the growth of bacteria and better protect your health.
- When not using this product for a long time, please keep the water tank dry.

FCC Statement

Any changes or modifications not expressly Approved by the party responsible for compliance could void the user's authority to operate the equipment. This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:
(1) This device may not cause harmful interference. (2) This device must accept any interference received, including interference that may cause undesired operation.
FCC Radiation Exposure Statement:
This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment.

WEEE Compliance

Correct Disposal of This Product (Waste Electrical & Electronic Equipment) (Applicable in countries with separate collection systems)
This marking on the product, accessories or literature indicates that the product and its electronic accessories should not be disposed of with other household waste at the end of their working life. To prevent possible harm to the environment or human health from uncontrolled waste disposal, please separate these items from other types of waste and recycle them responsibly to promote the sustainable reuse of material resources. Household users should contact their retailer, where they purchased this product, or their local government office, for details of where and how they can take these items for environmentally safe recycling. Business users should contact their supplier and check the terms and conditions of the purchase contract. This product and its electronic accessories should not be mixed with other commercial wastes for disposal.

FC FDA CE RoHS UK

Made in China

FCC Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement. The device can be used in portable exposure condition without restriction.