


User Manual Recupe sensor

1. Connect your sensor to your phone and the Recupe app (download from the App store).

With Recupe app open, pick up your Sensor. The sensor is motion activated and turns on when you pick it up and move it. When it does not turn on, give it a light shake in your hand. The sensor light should start flashing **green**.



2. After a couple seconds, you will see the words “Sensor Found” at the bottom of your screen. Your sensor is now connected to the Recupe app.



3. Hover the sensor over adhesive and it will stick to adhesive with the magnet.

4. When the Recupe sensor does not connect and light remains off, the battery is empty. Do not replace the battery yourselves but return it to Plety ; Use the grey bag which was provided when you received the sensor.

5. For product communication please include the serial number. The product serial number can be found on the label



6. FCC ID: 2BFOR-RECUPE

FCC Warning

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

NOTE: This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

The device has been evaluated to meet general RF exposure requirement.

The device can be used in portable exposure condition without restriction.