  I Model: K9244 V1.0 EN DE	<b>Table of contents</b> EN .....1-17 DE .....18-28	<b>1.Intended Use</b> <i>The device is not a substitute for medical consultation and treatment. Consult your doctor first if you are experiencing any pain or are suffering from an illness. The device is used for temporary relief of pain associated with sore and aching muscles in the waist or back, due to strain from exercise or normal household and work activities.</i> <b>2.Product Specifications</b> <table border="1"><thead><tr><th>Attribute</th><th>Specification</th></tr></thead><tbody><tr><td>Product Size</td><td>120x20cm / 47.24x7.87inch</td></tr><tr><td>Operating Power</td><td>Maximum 14.8W</td></tr><tr><td>Input Voltage</td><td>5V 2A</td></tr><tr><td>Battery</td><td>Rechargeable lithium-ion battery, 7.4V DC, 2000mAh, 14.8Wh</td></tr><tr><td>BLE Version</td><td>BLE 5.3</td></tr><tr><td>BLE Name</td><td>Comfytemp K9244</td></tr><tr><td>Charging Time</td><td>Approximately 3 hours (using 5V 2A adapter)</td></tr><tr><td>Endurance Time</td><td>Ranges from 90 - 480 mins depending on usage</td></tr></tbody></table>	Attribute	Specification	Product Size	120x20cm / 47.24x7.87inch	Operating Power	Maximum 14.8W	Input Voltage	5V 2A	Battery	Rechargeable lithium-ion battery, 7.4V DC, 2000mAh, 14.8Wh	BLE Version	BLE 5.3	BLE Name	Comfytemp K9244	Charging Time	Approximately 3 hours (using 5V 2A adapter)	Endurance Time	Ranges from 90 - 480 mins depending on usage	<b>3.Product Function</b> <p><b>A. Two Ways to Control:</b> Use the BLE App for convenient control or opt for onboard button control.</p> <p><b>B. Heating Function:</b> Choose from 9 levels.</p> <p><b>C. Vibration Function:</b> Choose from 5 vibration modes with 1~3 intensity level.</p> <p><b>D. Timer Function:</b> (1)Onboard Button Control: Default 30 minutes. (2)BLE APP Control: Choose a duration from 10 to 60 minutes.</p> <p><b>E. Battery Reminder:</b>Get notified when battery is low.</p> <b>4.Device Description</b> <ul style="list-style-type: none"><li>Power on/off</li><li>Heating button</li><li>Vibration button</li><li>Vibration Mode</li><li>BLE successful status</li><li>Heating Level</li><li>Low battery warning or charging status</li></ul>	<b>5. Use</b> <p><b>5.1 Wearing</b> When you receive the product, please charge it before use using a 5V 2A adapter for 2 hours. This ensures it is ready for proper functioning.</p> <ul style="list-style-type: none"><li>Tighten the outer strap.</li><li>Tighten the inner strap.</li><li>Tighten the outer strap.</li><li>Power on can be controlled with the button or the app.</li></ul> <p>To prevent potential skin burns due to individual sensitivity differences. <b>Please use over thin clothing, start at the lowest temperature setting for 20 minutes, and check skin condition after use.</b></p>	<b>5.2 Onboard Button Control</b> <p>(1) Step 1: Briefly press  to power on. (2) Step 2: Briefly press  to select the heating level. (3) Step 3: Briefly press  to choose the vibration mode. (4) Step 4: Long press  to select vibration intensity.</p> <p><b>5.3 BLE APP Control</b><p>(1) Scan the QR code to download the App or APK file. Or, search for "Comfytemp" in the app store to download it.</p><p>(2) Create an account.</p><p>(3) Briefly press to turn on the device, the BLE is ready to be connected. When BLE is connected, the icon  is on.</p></p>	<b>(4) Look for the device's BLE connection.</b> <p>(5) Begin using the products and associated apps.</p>	<b>6. Interaction (LED, Button, UI)</b> <table border="1"><thead><tr><th>Button</th><th>Operation</th><th>Function</th></tr></thead><tbody><tr><td>Heating Button</td><td>Briefly press</td><td>Choose from 3 levels (low / mid / high)</td></tr><tr><td>Power Button</td><td>Briefly press</td><td>Power on / off</td></tr><tr><td>Vibration Button</td><td>Briefly press</td><td>Choose from 3 vibration mode</td></tr><tr><td>Vibration Button</td><td>Long press</td><td>Choose from 3 vibration level</td></tr></tbody></table> <p>Note: Even when connected to the APP, selecting functions through the buttons only allows for 9 levels of heating and 3 vibration modes. For more functions, use the APP control.</p>	Button	Operation	Function	Heating Button	Briefly press	Choose from 3 levels (low / mid / high)	Power Button	Briefly press	Power on / off	Vibration Button	Briefly press	Choose from 3 vibration mode	Vibration Button	Long press	Choose from 3 vibration level	<b>6.2 UI</b> <table border="1"><thead><tr><th>UI</th><th>Functions</th></tr></thead><tbody><tr><td></td><td>BLE prompt</td></tr><tr><td></td><td>Temperature level display</td></tr><tr><td></td><td>Vibration mode selection</td></tr><tr><td></td><td>White: 1 vibration level Yellow: 2 vibration level Red: 3 vibration level</td></tr><tr><td></td><td>Fault alert F1 (for details, see section 10)</td></tr><tr><td></td><td>Low battery alert</td></tr><tr><td></td><td>Charging</td></tr><tr><td></td><td>Fully charged</td></tr></tbody></table>	UI	Functions		BLE prompt		Temperature level display		Vibration mode selection		White: 1 vibration level Yellow: 2 vibration level Red: 3 vibration level		Fault alert F1 (for details, see section 10)		Low battery alert		Charging		Fully charged
Attribute	Specification																																																										
Product Size	120x20cm / 47.24x7.87inch																																																										
Operating Power	Maximum 14.8W																																																										
Input Voltage	5V 2A																																																										
Battery	Rechargeable lithium-ion battery, 7.4V DC, 2000mAh, 14.8Wh																																																										
BLE Version	BLE 5.3																																																										
BLE Name	Comfytemp K9244																																																										
Charging Time	Approximately 3 hours (using 5V 2A adapter)																																																										
Endurance Time	Ranges from 90 - 480 mins depending on usage																																																										
Button	Operation	Function																																																									
Heating Button	Briefly press	Choose from 3 levels (low / mid / high)																																																									
Power Button	Briefly press	Power on / off																																																									
Vibration Button	Briefly press	Choose from 3 vibration mode																																																									
Vibration Button	Long press	Choose from 3 vibration level																																																									
UI	Functions																																																										
	BLE prompt																																																										
	Temperature level display																																																										
	Vibration mode selection																																																										
	White: 1 vibration level Yellow: 2 vibration level Red: 3 vibration level																																																										
	Fault alert F1 (for details, see section 10)																																																										
	Low battery alert																																																										
	Charging																																																										
	Fully charged																																																										

<b>7. Cleaning Instructions</b> <p>(1) Always keep the item away from water or other liquids, as this may cause damage. To clean the item, use a dry, lint-free cloth. Supported to use wet towel to wipe the heating pad surface. Not supported to wash the item. (2) Do not dry-clean this heating pad. (3) Do not use any chemical or abrasive cleaning agents. Cleaning solvents may have a deteriorating effect on the insulation of the heating element.</p>	<b>8. Storage</b> <p><b>When not in use, store as follows:</b></p> <ul style="list-style-type: none"><li>Allow the appliance to cool before folding when storing.</li><li>Avoid creasing the appliance by placing items on top of it during storage.</li><li>Regularly inspect the appliance for signs of wear or damage. If such signs are present, or if the appliance has been misused or is not functioning properly, return it to the supplier before attempting to use it again.</li><li>Do not fold the part with the heating sheet in half when storing, as it may damage the internal parts. Store in the original packaging, in a cool, dark, dry place.</li></ul>	<b>9.Danger</b> <p><b>IMPORTANT SAFETY INSTRUCTIONS</b> <b>DO NOT DESTROY</b></p> <p><b>WARNING – TO REDUCE RISKS OF BURNS, ELECTRIC SHOCK AND FIRE, THIS PRODUCT MUST BE USED IN ACCORDANCE WITH THE FOLLOWING INSTRUCTIONS:</b></p> <ol style="list-style-type: none"><li>READ ALL INSTRUCTIONS.</li><li>DO NOT USE WHILE SLEEPING.</li><li>DO NOT USE ON AN INFANT OR ON AN ANIMAL.</li><li>THIS PAD IS NOT TO BE USED BY OR ON AN INVALID, SLEEPING OR UNCONSCIOUS PERSON, OR A PERSON WITH POOR BLOOD CIRCULATION UNLESS CAREFULLY ATTENDED.</li><li>THIS APPLIANCE MUST NOT BE USED BY PERSONS INSENSITIVE TO HEAT AND OTHER VERY VULNERABLE PERSONS WHO ARE UNABLE TO REACT TO OVERHEATING.</li><li>THIS APPLIANCE CAN BE USED BY CHILDREN AGED FROM 8 YEARS AND ABOVE AND PERSONS WITH REDUCED PHYSICAL, SENSORY OR MENTAL CAPABILITIES OR LACK OF EXPERIENCE AND KNOWLEDGE IF THEY HAVE BEEN GIVEN SUPERVISION OR INSTRUCTION CONCERNING USE OF THE APPLIANCE IN A SAFE WAY AND UNDERSTAND THE HAZARDS INVOLVED.</li><li>CLEANING AND USER MAINTENANCE SHALL NOT BE MADE BY CHILDREN WITHOUT SUPERVISION. CHILDREN UNDER THE AGE OF THREE ARE NOT TO USE THIS APPLIANCE DUE TO THEIR INABILITY TO REACT TO OVERHEATING.</li></ol>	<p>8. BURNS CAN OCCUR REGARDLESS OF CONTROL SETTING, CHECK SKIN UNDER PAD FREQUENTLY. 9. DO NOT USE PINS OR OTHER METALLIC MEANS TO FASTEN THIS PAD IN PLACE. 10. DO NOT SIT ON, OR CRUSH PAD. AVOID SHARP FOLDS. 11. UNPLUG WHEN NOT IN USE. NEVER LEAVE APPLIANCE UNATTENDED ESPECIALLY IF CHILDREN ARE PRESENT. 12. CAREFULLY EXAMINE PAD, COVER AND CONTROLLER BEFORE EACH USE. DISCARD THE PAD IF ANY SIGN OF DETERIORATION, WEAR OR DAMAGE. DO NOT TAMPER WITH THIS PAD IN ANY WAY. THERE ARE NO USER SERVICEABLE PARTS. IF PAD DOES NOT FUNCTION SATISFACTORILY, SEE WARRANTY FOR CUSTOMER CONTACT INFORMATION. 13. USE THIS PAD ONLY AT RATED VOLTAGE, AS MARKED ON THE PRODUCT. DO NOT USE WITH GENERATORS, POWER CONVERTERS, OR INVERTERS. 14. DO NOT USE THIS PAD WITH LINIMENT, SALVE OR ORNIMENT PREPARATIONS. SKIN BURNS COULD RESULT. 15. THIS APPLIANCE IS NOT INTENDED FOR MEDICAL USE IN HOSPITALS. 16. DO NOT USE IF WET (NOT REQUIRED FOR MOST WET-PROOF OR CLASS III APPLIANCES). 17. SAVE THESE INSTRUCTIONS. 18. PLEASE CHECK THIS HEATING PAD AND CONNECTION CABLE FREQUENTLY FOR SIGNS OF WEAR AND TEAR OR DAMAGE. IF DETECTED, DISCONTINUE USE AND CONTACT CUSTOMER SERVICE. 19. WHILE CHARGING, DO NOT COVER THE POWER BANK WITH COVERS, CUSHIONS OR SIMILAR ITEMS. 20. DO NOT CHARGE THE POWER WHILE HEATING PAD IS IN USE. 21. RISK OF EXPLOSION! DO NOT PUT POWER INTO FIRE OR IT MAY EXPLODE.</p>	<b>10. Trouble Shooting</b> <p><b>Problem01:</b> The device does not switch when the ON/OFF button is pressed. <b>How to proceed:</b> The device runs out of battery. Please charge with a 5V 2A adapter for 3 hours before using again. When the LED stays green, it means it is fully charged.</p> <p><b>Problem02:</b> The device shut down after working for a while. <b>How to proceed:</b> (1) The current voltage of the battery is too low. Please make sure to charge with a 5V2A adapter for 3 hours until fully charged before use. (2) When BLE is not connected, default timer is 30 minutes. If longer timer settings are needed, please use the app.</p> <p><b>Problem03:</b> The device is not hot enough. <b>How to proceed:</b> (1) It is recommended to use it continuously for more than 10-20 mins before judging if the temperature meets the needs. (2) Low battery levels affect temperature performance. If low battery is detected, please recharge promptly. (3) In order to prevent low temperature burns, we have limited the highest temperature.</p> <p><b>Problem04:</b> Skin in the treatment area turns red. <b>How to proceed:</b> If the redness disappears quickly, there is no risk – this is caused by the locally compressed, increased blood flow. Here is some suggestion: • Moderate heat: Choose an appropriate heat level to avoid skin irritation. • Use protection: Place a thin cloth between skin and heat pad for added comfort. • Time and frequency: Limit use to 20-30 minutes per session, with breaks for skin recovery. • Avoid during inflammation: Don't use on inflamed skin to prevent worsening symptoms.</p> <p><b>Problem05:</b> Can't connect to BLE. <b>How to proceed:</b> (1) First, make sure your phone's BLE is turned on. (2) Follow the instructions on Section 5.3 of User Manual for BLE connection. (3) Try disabling the BLE device and set it up again. (4) For safety, you can't use the device while charging the product. (5) To find the problem, check if your phone can connect to other BLE devices at home or try a different phone to connect to this one.</p> <p><b>Problem06:</b> The phone can detect BLE, but cannot connect to it. <b>How to proceed:</b> For Android system, please go to Settings-App &amp; permissions-App Info-Top right corner (three dots)-Show system apps-BLE-Storage- Clear data, then attempt to reconnect.</p>	<b>Warranty</b> <p>24 Hrs Quick Customer Support 30 Days Free Return 365 Days Warranty</p> <p>If you have questions or concerns, don't hesitate to contact us directly at: US Customer Services Hotline <b>1(833)990-0618</b> Monday to Friday, 9 AM to 5 PM EST Email: <b>help@comfytemp.com</b></p>	<b>1.Vorgesehene Verwendung</b> <p>Das Gerät ersetzt keine ärztliche Beratung und Behandlung. Konsultieren Sie zuerst Ihren Arzt, wenn Sie Schmerzen verspüren oder an einer Krankheit leiden. Das Gerät dient zur vorübergehenden Linderung von Schmerzen, die mit verspannten und schmerzenden Muskeln in der Taille oder im Rücken aufgrund von Belastung durch Bewegung oder normale Haushalts- und Arbeitsaktivitäten verbunden sind.</p>
--	--	--	--	--	---	--

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**UK REP** Company: YH Consulting Limited  
Address: C/O YH Consulting Limited Office 147, Canton House, London Road, Slanes upon Thames, Surrey, TW18 4AX, London  
E-mail: H27HJK@gmail.com  
Telephone: +44 07514 677868

Manufacturer: Shenzhen Nua Health Technology Co., Ltd  
Address: 11th Floor, Zhenshen Building, Yousong Community, Longhua Street, Longhua District, Shenzhen, Guangdong, China, 518110.

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**FCC Warning:**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

**Caution:** Any changes or modifications to this device not explicitly approved by manufacturer could void your authority to operate this equipment.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

This equipment complies with FCC radiation exposure limits set forth for an uncontrolled environment. This equipment should be installed and operated with minimum distance 0cm between the radiator and your body.

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**EU REP** Company: E-CrossStu GmbH