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To charge the product, first open the charger port cover at the bottom of the handle. Connect the charger cable that comes along with the product to a USB charger adapter in line with national standards. Place the electronic toothbrush in a dry place for charging. It takes approximately four to five hours to fully charge the unit Charger port cover AS03 Portable Sanitizer Battery Status and Charging On/Off Button/Mode Button Explanation of Brushing Modes: How to Install a Special usmile Brush Head Operating Instructions for the Electroni-This product and the portable sanitizer are charged prior to Toothbrush: Please use the product together with a special usmile brush head. Replace the special usmile brush head after three months of use or when the bristles start to fall out or split so as to prevent discomfort to your gums and teeth. (the duration cannot be adjusted in the app). Before the timed shutdown, the frequency will change three times as a reminder. Special brush head x 2 Charger cable x 1 Charger port **Explanation of Customized Forces:**  Steplessly adjust the customized forces of all modes in the brush function settings in the app. Handle of the usmile F1 Smart ElectronicTooth- Switch between forces: Under the power-on mode, short press the force button to switch between forces Operating Instructions of the Portable Sanitizer: **User Manual** 31 mm x 31 mm x 253 mm (including the brush head) Red: Low battery. Please charge. Motor shaft place the charger cable in water or a humid environment. riging, be sure the charger cable and the charging of the product are dry and clean. Stains or impurities can affect the safety of charging. riging, please detach the charger cable from the product off the power. Open the charger port cover at the bottom **Future One Smart Sonic** Force indicator Adjust the Brushing Mode and Start to Electric Toothbrush Brush Your Teeth On/Off button in-one charger cable, which can damage the Force button Red: Low battery. Please charge. Mode button The first usmile indicator is on indicator is on usmile Clean mode
White mode
Soft mode
Massage mode
Smart mode Smart mode-----Power/Bluetooth indicator-----Power/Bluetooth indicator Power indicator light Charging interface The third usmile indicator is on All three indicators are on usmile Be sure you have covered the charger port after charging before you start to use the electronic toothbrush. \* For the smart mode, you can select the "Polish Mode、Gum ca Mode 、Cigarette stain Mode and Tongue Mode"in the app . frequency will change once to remind you of shifting to another section of your mouth. Before the timed shutdown, the frequency Model: Future One Schematic Drawings of Components 5. Begin to brush each surface (see the figures below) Travel with the Product Correct Tooth Brushing Method Cleaning and Storage 4. Hold the handle and press the shaft down to push the movement component out of the handle. Install the usmile Oral Care No other phones were linked with the current toothbrush. Fail to insert the brush head into the motor shaft motor shaft motor shaft ı have already downloaded the app, be sure it is the latest version.

Follow app instructions to add a device. It is suggested that when brushing the outer and inner surfaces, you hold the brush at an angle toward the gum line, as shown in the figure below, to better clean your teeth.

Keep the electronic toothbrush in a stable position to prevent it from tumbling or falling.

The postable satisface is not watercross? Place it in a well-ventilated.

The postable satisface is not watercross? Place it in a well-ventilated. Please be sure your toothbrush has battery life. Place the toothbrush close to your mobile phone. 2.Go to the home screen of the app. The app will pop up a b requesting connection. Follow instructions to add the device. Not a usmile brush Replace it with a usmile brush head Treat . 8. If there is no such pop-up box after going to the home screen of app, click "+" in the upper right corner to manually add the devic instructed.
>>When the power/Bluetooth indicator turns blue, it indicates the toothbrush has been connected to the app. Environmental Protection The brush head vibrates weakly 4.Upon connection, you can adjust the mode and force of the toothbrushin the app. 1.Why does the Android app request location permissions? For Android Version 6.0 and above, the system requests users to grant location permissions before using the app.

1.Why does the Android app request location permissions? For Android Version 6.0 and above, the system requests users to grant location permissions before using the app. **Faults and Processing Methods** Charging

The working time after charging is too short

The charging time is too short

The charging time after charging time at the constant the constant that the constant the constant that the product indoors at 15°C-35°C

The charging time is too short

The charging time is too short

Charge the product constant the device and the product indoors at 15°C-35°C

The charging time is too short

The working time perature is too low product indoors at 15°C-35°C

The charging time is too short

Charge the product indoors at 15°C-35°C

The device has been evaluated to meet general RF exposing time requirement. The device can be used in portable exposu condition without restriction. Fault Possible Reasons Solutions **-**() 2.Possible reasons for the failure to find a Bluetooth device include: Fail to turn it on by pressing the On/Off button

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