## **Pendant FD User Guide**



Medical Alert Systems providing

Independence and peace of mind.


## **Pendant FD User Guide**

### 


### **How Pendant FD works**

#### **Use Case 1 : EMERGENCY BUTTON press**


1. Press the emergency button of pendant FD.



3. Emergency operators establishes two-ways communication and assesses the situation.



2. Your information and location are sent to the Emergency Operators.



4. Emergency operator contacts a neighbor, family member or emergency services based on specific need.



### **Pendant FD User Guide**

#### **Use Case 2: Fall Detection**

When you experience a fall in the shower or car park anywhere you are out of reach of a cell phone, pendant FD can send the signal for help in faster.



#### Fall detect with and without cancellation

- 1. If you experience a fall, pendant FD waits twenty seconds to confirm that the user did not get back up. Pendant FD broadcasts Fall Detect RF signal. If no cancellation is received after twenty seconds, pendant FD will establish full alarm. The alarm is sent to Emergency Operators. Emergency operator will assist based on specific need.
- 2. If you experience a fall, pendant FD waits twenty seconds to confirm that the user did not get back up. Pendant FD broadcasts Fall Detect RF signal. If you presses and holds the Pendant for five seconds, the LED will be turn off and the alarm will be cancelled.

# Freedom to live life your way

### **Federal Communication Commission Interference Statement**

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

FCC Caution: Any changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate this equipment.

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.