Silk'nTM
Beautiful Technology

Lipo

H5001
H5002

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
Table of contents

1. Warnings & contraindications.....	2
2. Safety	8
3. Get to know Silk'n Lipo.....	8
4. Device controls & settings.....	10
5. Treatment plan.....	15
6. What to expect with Silk'n Lipo	16
7. Get started!.....	17
8. Cleaning & maintenance	19
9. Parts and consumables	20
10. Troubleshooting	21
11. Customer service	21
12. Warranty	21
13. Disposal	22
14. Specifications	23
15. Signs & symbols	25
16. Labeling	26
17. Compliance Statements	27

1. Warnings & contraindications

EN

Please read all warnings and contraindications before use.

-  Do not use the device while bathing or taking a shower.
-  Always ensure that your hands are dry before connecting the stand/charging cradle to the electricity supply.
-  If the stand/charging cradle malfunctions, unplug immediately.
-  Dry the electrode pad/s before storage.
-  Only authorized personnel are permitted to perform repairs. Do not attempt to modify, open or repair the device. You may expose yourself to dangerous electrical components and may cause serious harm to yourself.
-  Do not use if the parts or accessories are not supplied by the manufacturer as defined in this manual, or if the device is damaged or does not appear to be working properly. In such cases contact Customer Service.

- ⚠ Changes or modifications to this device not expressly approved by the party responsible for compliance could void the user's authority to operate the device.
- ⚠ Do not clean any of the parts in the dishwasher.
- ⚠ Do not dispose with household waste! The law requires that you dispose of this electrical device (including all accessories and parts) at a designated recycling collection point for electrical devices. For details see section "13. Disposal").
- ⚠ The battery must be removed from the appliance before it is scrapped.
- ⚠ The appliance must be disconnected from the supply mains when removing the battery.
- ⚠ The battery is to be disposed of safely.
- ⚠ Attention should be drawn to the environmental aspects of battery disposal. Don't throw used batteries in a dustbin. Please contact your retailer in order to protect the environment.

Battery Removal

EN

To remove the battery, the device must be broken apart. The device will no longer function after this procedure, nor can it be repaired. Any warranty will be void when the device housing is broken and the battery is removed. We recommend you wear safety protective goggles and gloves.

Silk'n is not liable for any physical or material damages as a result of performing this procedure.

To remove the battery place the hand unit on a hard, strong surface. Use a hammer to hit the device on the seam and continue to do so along the whole seam on both sides of the device.

Once the seam is sufficiently broken, use pliers or a screwdriver to take the housing apart. Unscrew the PCB (printed circuit board) and remove the battery by disconnecting the plug from the PCB.

Contraindications

- ⚠ Silk'n Lipo must not be used by children, adults with limited physical, sensory or psychological capacities, or those lacking experience or knowledge in how to use the device. Keep the device out of the reach of children.
- ⚠ Use a lower EMS energy level if you experience any discomfort.
- ⚠ Stop use immediately if your skin is very red, or if it blisters or burns.
- ⚠ Use on clean, dry, healthy skin.

Do not use if you:

- ⚠ have heart disease and you have not consulted your physician.
- ⚠ are on immuno-suppressant drugs such as corticosteroids without consulting your physician, or if you have an immuno deficiency disorder.
- ⚠ have any active type of cancer.
- ⚠ have a lymphatic disorder, thyroid disorder, diabetes or liver disorder, or any other metabolic condition.

- ⚠ have a pacemaker, internal defibrillator, or any other active electrical implant anywhere in the body.
- ⚠ have applied topical solutions containing metal ions on the treatment area.
- ⚠ have epilepsy.
- ⚠ have serious arterial circulatory problems in the lower limbs.
- ⚠ have an abdominal or inguinal hernia.
- ⚠ are pregnant or nursing.

Do not use on these parts of the body:

- ⚠ Near the heart. The electrodes must not be placed on any part of the front ribcage (where the ribs and breastbone are located), and specifically not on the two large pectorals. This can increase the risk of ventricular fibrillation and induce cardiac arrest.
- ⚠ The back ribcage.
- ⚠ The head or neck.
- ⚠ On the genitals.
- ⚠ On or around the eyes.

- ⚠ Areas that have metabolic arthritis or chronic pain.
- ⚠ Areas that are swollen, infected, inflamed areas, broken or irritated skin, or skin eruptions (e.g. skin rashes, phlebitis, thrombophlebitis, varicose veins, or post-operative scars where the healing process could be affected etc.).
- ⚠ Parts of the body that have metal staples, implants or pins; or if you are near metal. Remove jewelry or other accessories from the treatment area during a session.

2. Safety

With Silk'n safety comes first. Silk'n Lipo has these safety features:

- **Proximity sensor:** The device works only when it is on, and the treatment unit/s are in good contact with your skin.
- **Auto shut-down:** Although you can stop a treatment whenever you wish, the device shuts down after a minute or so, when there is no contact with the skin (the *operating light* will be constant green). It also shuts down automatically after 15 minutes of continuous use.

3. Get to know Silk'n Lipo

Silk'n Lipo is a system to help you reduce fat in specific areas, exercise your muscles and tone specific parts of your body without physical effort, and in the comfort of your home. By placing the treatment unit/s on target areas such as the upper arms, abdomen, buttocks or thighs, fat can be broken down and removed specifically from that area. Silk'n Lipo is not intended for overall body fat reduction, and cannot replace diet or exercise, but unlike diet or exercise, Silk'n Lipo can help you shape individual areas.

Silk'n Lipo is a powerful electrical device and should be used with special attention to safety. **Before use, it is important that you read and understand the instructions on how to use the device including after-use procedures, and follow them strictly.** We recommend you re-familiarize yourself with the instructions in this manual before each treatment.

3.1 Package contents

The Silk'n Lipo pack includes:

- 1 or more Treatment units *
- 2 Connectors
- 1 to 3 Body belts *
- 1 or more Electrode pads *
- Stand/charging cradle + adapter
- Quick guide
- User manual
- Warranty card

*Varies according to the option purchased.

3.2 Device description

Silk'n Lipo is a device for home use that utilizes Low Level Laser Technology (LLLT) and Electrical Muscle Stimulation (EMS). Together these energies help increase muscle tone and reduce fatty acids and deposits to enhance body toning and reshaping.

Silk'n Lipo includes a selection of belts suited to the part of the body you wish to treat. Each treatment unit contains 8 laser diodes and electrodes that activate the EMS. The treatment unit/s are attached to the belt and then placed on the part of the body you wish to tone and reshape. Silk'n Lipo uses rechargeable batteries so the treatment unit/s can be easily charged before each use.

Refer to the diagram at the front of the manual to see the location of each feature.

- | | |
|---------------------|------------------------------|
| 1. Treatment unit | 5. Body belt |
| 2. Control button | 6. Electrode pads |
| 3. Indicator lights | 7. Charging indicator lights |
| 4. Connectors | 8. Stand / charging cradle |

3.3 Intended use

Silk'n Lipo is a non-invasive aesthetic treatment to reduce the circumference of the waist, arms, thighs and buttocks. The device is also intended to improve muscle tone and firmness, and strengthen muscles in these areas.


The device is not intended for use in any therapy, or for the treatment of any medical conditions or diseases. Do not use for any purpose other than as described in this manual.

4. Device controls & settings

4.1 Connecting the different parts


Connecting Silk'n Lipo is easy. First connect the *treatment unit/s* to the *belt*, and then wear the *belt* on the part of the body you wish to treat.

1. Remove the different parts from the box.
2. Place the *treatment unit/s* on a flat surface. If you have more than one *treatment unit*, connect them together.
3. Take the female *connector* and slot it into the left side of the *treatment unit*. Then take the male *connector* and connect it to the right side.


4. Select the belt for the part of the body you are treating, and lay it on a flat surface with the light gray side facing you.

5. Next, fold each side of the belt inwards so that you see the black Velcro pads.

6. Slide each end of the belt into the *connectors* and then fold backwards, so that the Velcro pad sticks to the belt. Your belt will now look like this.


7. Of course, to use Silk'n Lipo, you will need to disconnect one side, place the *belt* on your body and then re-connect the side you disconnected!

Note that you can only use the device when:

- the *treatment unit/s* are charged.
- you have placed a wet *electrode pad* on each *treatment unit* (see section "4.4 Wetting the electrode pads").
- you have read this manual, and understand how to use the device.

4.2 Controls and settings

Battery / Error Operating Not in use

Each *treatment unit* has a control panel. The control panel includes:

- a control button (#2)
- 2 EMS energy level buttons (left and right of the control button)
- 2 indicator lights - left and middle (#3).

Control	Action	Function
Control button 	Long press	• Turns the <i>treatment unit</i> on (or off). When you turn it on, you will feel a short vibration. • The green <i>operating light</i> (middle) will blink continuously.
+ Button 	Short press	• Each press increases the EMS level. As you increase levels, the <i>operating light</i> blinks faster.
	Long press	• The same effect as a short press.
- Button 	Short press	• Each press decreases the EMS level.
	Long press	• Decreases the EMS level to zero.
Battery / Error light (left)	Orange light blinks continuously	• Battery low.
	Constant orange	• Device error.

Control	Action	Function
Operating light (middle)	Green light blinks continuously	• The device is on, and has good contact with the skin. • Blinks faster as the EMS level increases.
	Constant green	• There is no contact with the skin. • The device will then shut down automatically after a minute or so.
Indicator light (right)	---	• Not in use.

4.3 Energy levels

Silk'n Lipo has multiple levels of Electrical Muscle Stimulation (EMS). This gives you excellent control of the intensity of muscle stimulation you wish to use for each part of the body you treat. What's more, you can change levels as often as you like during a treatment.


- **To increase the EMS level:** Short press the *+* button as many times as you need to set to the EMS level that is comfortable for you. You will see that *operating light* blinks more quickly as you increase the EMS level.
- **To reduce the EMS level:** Short press the *-* button as many times as you need (refer to section "4.2 Controls and settings").
- **To go back to the lowest EMS level:** Long press the *-* button.

i You can only change EMS levels when the *treatment unit/s* are in good contact with your skin. When there is no contact with the skin, the EMS drops to zero.

4.4 Wetting the electrode pads

The *electrode pads* are required for the EMS to work, and they must be wet before you start your treatment.

1. Wet each *electrode pad* with a wet cloth or with wet fingers.
2. Then place the wet *electrode pad* onto the underside of each *treatment unit*. Make sure that it fits firmly around the center rectangle. You might need to manipulate it slightly, so that the *electrode pad* lies flat on top of the *treatment unit*.

i Only wet the *electrode pad*, not the whole *treatment unit*.

4.5 Device charging

Charge the *treatment unit/s* before first use. The initial charge should take around 6 hours. Subsequent charges should take around 3 hours. We recommend you re-charge the device before each treatment.


1. Remove the **electrode pads** before placing the **treatment units** in the **stand/charging cradle**.
2. Connect the empty *stand/charging cradle* to the electricity supply. The 2 *charging indicator lights* will be constant green.
3. Place the *treatment unit/s* in the cradle so that they sit in the designated grooves.
 - When the *treatment unit* is charging, the orange *battery/error indicator light* on the *treatment unit* will flash on and off.
 - When fully charged the *battery/error indicator light* (left) will go off, and the green *operating light* (middle) will be constant green.

4. Once charged, disconnect the *stand/charging cradle* from the electricity supply.

5. Treatment plan

EN

- ⚠ For use on clean, healthy skin (see section "1. Warnings & contraindications").
- Treat the target area for 15 minutes per treatment session.
- For larger parts of the body, like the stomach or buttocks, you may wish to move the *treatment unit/s* to an adjacent part of the target area, and use for another 15 minutes. You may keep doing this until you have covered the target area.
- Try to avoid gaps and overlaps.
- Never use Silk'n Lipo on exactly the same spot during a treatment session.
- For daily use, but no more than 3 times a week on each specific target area. Continue until you get the results you want.

6. What to expect with Silk'n Lipo

If you follow the instructions fully and pay attention to **all warnings and contraindications**, complications and side effects are uncommon. However every energy-based device, including those designed for home use, involves some degree of risk. Therefore it is important that you understand and accept the risks and complications that may occur. Be sure to review section "16. Labeling".

Skin redness

You may find some skin redness in the treated area which should disappear after the treatment.

Muscle pain

As with all exercise, some muscle soreness may occur after a treatment but it should go away after a day or two. If this happens, use a low EMS level for a few sessions until the soreness has gone.

Skin irritation

Although very rare, it is possible to have a prolonged reddening of the skin, or the appearance of acne. Stop use immediately, and only treat that specific area when the skin irritation has gone.

Headache

Equally rare, is the possibility that you might feel light-headed when you start using the device. If this happens, sit down during the treatment, or stop using the device until you feel OK.

Allergic reaction

If you have an allergic reaction to the device or its accessories, stop using the device.

 Stop use immediately if any of these effects are more severe.

 See your doctor if the effects do not go away within 2 to 3 days.

7. Get started!

You're on your way to a trim and toned body! Treating with Silk'n Lipo is easy. Just follow these steps for each part of the body you wish to treat.

7.1 Treatment steps

1. Use on clean, healthy skin (refer to section "16. Labeling").
2. Check that the batteries are well charged before use. If the *battery/error light* is orange and blinks slowly, this indicates that the battery is low.
3. Before each treatment, wet the *electrode pad/s* with a wet cloth or with wet fingers, and place them on each *treatment unit* (see section "4.4 Wetting the electrode pads"). Next, connect one or more *treatment unit/s* to the *belt* (see "4.1 Connecting the different parts").
4. Disconnect one of the *connectors* on one side, and then place the *belt* around the area you wish to treat. Then reconnect the side you just disconnected.
5. Adjust the tightness of the *belt* so that the *treatment unit/s* are in good contact with your skin, and that the *belt* is firm and well secured to your body.

 If the *belt* is not tight against your skin, the EMS energy may disconnect. If that happens, press the *+ button* to re-select the EMS energy level required.

6. **Switch on.** With the controls in view, press the *control button* on each *treatment unit*. The green and orange *indicator lights* will go on, and then go off after a second or so and the *treatment unit/s* will vibrate once. When the *treatment unit/s* are in good contact with the skin, the green *operating light* will blink slowly, and the laser energy will start to work.

7. Changing the EMS level:

- **To increase:** Short press the *+ button* as many times as you need to set to the EMS level that is comfortable for you. The *operating light* blinks more quickly as you increase the EMS level.

 If you long press the *+ button*, the EMS level will still only increase by one increment.

- **To reduce:** Short press the *- button* as many times as you need (refer to section “4.2 Controls and settings”).

 Long pressing the *- button*, turns off the EMS energy.

8. **Treat each area for 15 minutes.** Keep the *treatment unit/s* in place, and the device on for up to 15 minutes. Every few minutes you will feel a vibration. At the end of 15 minutes, you will feel 3 vibrations, no more energy will be emitted and the device will shut down automatically.

 To switch off the device **before 15 minutes**, long press the *control button*. You will feel 2 vibrations, then the device will shut down.

9. If you want to treat another area (next to the area you have just treated or another part of your body), wet the *electrode pad/s*, and move the *treatment unit/s* to the next location.

 If you move the device to another part of your body and the *treatment unit/s* are not in good contact with the skin, the EMS level will revert to zero. If this happens, make sure the *treatment pad/s* are in good contact with the skin, then short press the *+ button* to increase the EMS level.

8. Cleaning & maintenance

8.1 Cleaning the treatment units & charging cradle

We recommend that you clean the *treatment unit/s* after each use, and the *stand/charging cradle* once a month.

Before cleaning:

- Remove the *treatment unit/s* from the *stand/charging cradle*.
- Disconnect the *stand/charging cradle* from the electricity supply.
- Do not use abrasive cleaners or chemicals to clean any parts.
- Clean the *treatment unit/s* or *stand/charging cradle* with a damp cloth using a mild detergent, then wipe it dry with a clean cloth.

Do not clean any of the parts in the dishwasher.

8.2 Cleaning the electrode pads

- After each use, wipe the *electrode pad/s* and dry with a cloth.
- Do not store them in a closed box or drawer until they are completely dry.

8.3 Replacing the electrode pads

- Do not use *electrode pad/s* if they get damaged. If that happens, replace with new ones.

8.4 Cleaning the belt

Hand wash only. Please follow the instructions on the belt.

External: Nylon. Middle
layer: SBR. Internal: Nylon

8.5 Storage

- After cleaning, we recommend that you return the *treatment unit/s* to the *stand/charging cradle*.
- Store the different parts in the original packaging in a clean, dry place to prevent the accumulation of dust.
- Make sure the *stand/charging cradle* is disconnected from the electricity supply.

9. Parts and consumables

You can purchase replacement treatment unit/s or belts through the retailer where you bought the device, or you can order directly from www.silk'n.com.

10. Troubleshooting

Only authorized personnel are permitted to perform repairs. Do not attempt to open or repair your device. You may exposure yourself to dangerous electrical components and cause serious harm to yourself. It will also void the warranty.

Problem	Checks
My device does not start.	Make sure that: • the <i>treatment unit/s</i> are charged. • you have pressed the <i>control button</i> to turn them on. • the <i>operating light</i> is blinking.
I can't feel the electro stimulation or the EMS level seems too low.	Check that: • the <i>electrode pad/s</i> are wet. • the <i>treatment unit/s</i> are in good contact with the skin and • the <i>belt</i> is fixed securely to your body.
I can't change the EMS level.	Check that the <i>treatment unit/s</i> are charged.
The indicator lights are on all the time.	There is a system error.

 If you have problems or queries, please contact Customer Service.

11. Customer service

See our website for more information: www.silkn.com. This manual is also available as a PDF download from: www.silkn.com/manuals.

Contact Customer Service to report unexpected operations or events, or any other problem with your device. Our contact details are:

Email: contact@silkn.com Telephone: 1-877-367-4556 (toll-free)

12. Warranty

EN

Refer to your Silk'n warranty card.

13. Disposal

If your device reaches its End of Life stage, be sure to follow your country's laws for the disposal of electrical items.

Dispose of the packaging separated into single type materials. Dispose of cardboard and carton as waste paper and film via the recyclable material collection service.

Do not dispose with household waste! The law requires that you dispose of this electrical device (including all accessories and parts) at a designated recycling collection point for electrical devices. This ensures that the device is recycled professionally, and prevents the release of harmful substances* into the environment. Make sure that the device is in an uncharged state before disposal.

*Products labelled with Cd = cadmium, Hg = mercury, Pb = lead.

14. Specifications

Model No.	H5001 (device); H5002 (charging cradle)
Technology	Low Level Laser Therapy (LLLT) Electrical Muscle Stimulation (EMS)
Charging cradle	Input 12Vdc, 1.5Amax Charging 2 modules simultaneously
Module powered by	Li-ion battery, 3.7V, 2600mA
Adapter models	PMA18A-120 (USA) KSA24A1200150HU (USA)
Adapter Rating Input: Output:	100-240 Vac, 50-60 Hz, 0.6A Max. 12 Vdc, 1.5A
LLLT Diode wavelength Diode power	658nm +/- 8nm 100mW (max)
EMS Pulse Amplitude	0-32mA into 1000 ohm load
Output Voltage	Adjustable 0-32V, Max output 32V peak to peak into 1000ohm load.
On Time	Adjustable, 20~2400 milliseconds
Off Time	Adjustable, 20~2400 milliseconds
Ramp Time	Adjustable, 22~800 milliseconds, 22 milliseconds/step. The "On" time will increase and decrease.
Pulse Rate	30~50 Hz
Pulse Width	350~450µs
Timer	15 min
Wave Form	Symmetrical Bi-phasic rectangular pulse
Tolerance	There may be a +/-5% tolerance of all settings and +/- 20% tolerance of output of intensity.

Package size	192 x 296 x 130 (L x l x H) [mm]
Weight	Module: 300 [gr] Charging cradle: 215 [gr]
Transport & storage between uses and storage condition	Temperature: minus 40 - 70 [°C] Relative humidity: 10% - 90% Atmospheric pressure: 500 - 1060 [hPa]
Operating conditions	Temperature: 5 - 40 [°C] Relative humidity: 15% - 90% Atmospheric pressure: 700 - 1060 [hPa]
Time from minimum storage temperature between uses to minimum operating temperature with ambient temperature of 20°C degrees	30 minutes
Time from maximum storage temperature between uses to maximum operating temperature with ambient temperature of 20°C degrees	30 minutes
Service life	3 years

15. Signs & symbols

EN

	Declaration of Conformity: Products labeled with this symbol meet all applicable provisions of the European Economic Area.
	CSA mark for US and Canada.
	Assembled from tested components.
IP22	The device is protected against the ingress of objects such as fingers.
	Degree of protection against electric shock. BF applied part.
	Keep dry.
	Do not dispose with household waste! The WEEE directive requires that you dispose of this electrical device (including all accessories and parts) at a designated recycling collection point for electrical devices.
	Follow the operating instructions.
	Class 1 laser product.

16. Labeling

EN

HomeSkinovations
Tavor Building, Shaar Yokneam,
POB 533, Yokneam 2069206, Israel

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IP22

CE, FCC, RoHS, Laser 1

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference and (2) This device must accept any interference received, including interference that may cause undesired operation.

Made in ISRAEL

HomeSkinovations
Tavor Building, Shaar Yokneam,
POB 533 2069206 Yokneam, Israel

MODEL No.: H5002

2017-08-9

12V --- 1.5A

CE, RoHS, Laser 1

Made in ISRAEL

17. Compliance Statements

The use of accessories other than those specified or sold by Home Skinovations Ltd as replacement parts may have the consequence of increasing the emissions or decreasing the immunity of the device. The device is intended for use in the electromagnetic environment specified in tables 1, 2, 4 and 6.

Table 1: Guidance and Manufacturer's Declaration: Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment Guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	Complies The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Not Applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2: Guidance and Manufacturer's Declaration: Electromagnetic Immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the devices should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines NA	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	2 kV 0, 90, 180 & 270 degree phase angles	2 kV 0, 90, 180 & 270 degree phase angles	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	100% drop, 0.5 periods (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) 100% dip, 1 period 30% dip, 25/30 periods 100% drop, 5 seconds	100% drop, 0.5 periods (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) 100% dip, 1 period 30% dip, 25/30 periods 100% drop, 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note <i>UT is the a.c. mains voltage prior to application of the test level.</i>			

Table 4: Guidance and Manufacturer's Declaration: Electromagnetic Immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.

Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	6 Vrms 150 kHz to 80 MHz	6 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance. $d = \sqrt{\frac{3S}{P_1 N P}}$
Radiated RF	10 V/m	10 V/m	$d = \sqrt{\frac{3S}{E_1 N P}}$ 80 MHz to 800 MHz
IEC 61000-4-3	80 MHz to 2,7 GHz	80 MHz to 2,7 GHz	$d = \sqrt{\frac{7}{E_1 N P}}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

A Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.

B Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Table 6: Recommended Separation Distances between Portable and Mobile RF Communications Equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2,5 GHz
	$d = \left[\frac{3,5}{V_1}\right]\sqrt{P}$	$d = \left[\frac{3,5}{E_1}\right]\sqrt{P}$	$d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.7	3.7	7.38
100	11.7	11.7	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

FCC Compliance

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- a) Reorient or relocate the receiving antenna.
- b) Increase the separation between the equipment and receiver.
- c) Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- d) Consult the dealer or an experienced radio/TV technician.

RF Radiation Exposure

This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment. For hand-held/body-worn operation, this equipment has been tested and meets the FCC RF exposure guidelines. This transmitter must not be co-located or operated in conjunction with any other antenna or transmitter. Use of other accessories may not ensure compliance with FCC RF guidelines.

Notice for Canada

This device complies with Industry Canada licence-exempt RSS standard(s). Operation is subject to two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference that may be received or that may cause undesired operation.

A distance of at least 0.5cm between the equipment's antenna and all persons should be maintained during the operation of the equipment.

© 2017 Home Skinovations Ltd.

All rights reserved.

EN

Legal Notice

Home Skinovations Ltd reserves the right to make changes to its products or specifications to improve performance, reliability, or manufacturability. Information furnished by Home Skinovations Ltd is believed to be accurate and reliable at the time of publication. However, Home Skinovations Ltd assumes no responsibility for its use. No license is granted by its implication or otherwise under any patent or patent rights of Home Skinovations Ltd.

No part of this document may be reproduced or transmitted in any form or by any means, electronic or mechanical, for any purpose, without the express written permission of Home Skinovations Ltd. Data is subject to change without notification.

Home Skinovations Ltd has patents and pending patent applications, trademarks, copyrights, or other intellectual property rights covering subject matter in this document. The furnishing of this document does not give you any license to these patents, trademarks, copyrights, or other intellectual property rights except as expressly provided in any written agreement from Home Skinovations Ltd. Specifications are subject to change without notice.

SILK'N and the SILK'N logo are registered trademarks of Home Skinovations Ltd, Tavor Building, Shaar Yokneam POB 533, Yokneam 2069206, ISRAEL www.silkn.com, contact@silkn.com.

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INSTRUCTIONS FOR USE


Table of contents

1. Warnings & contraindications.....	2
2. Safety	8
3. Get to know Silk'n Model H5003	8
4. Device controls & settings.....	10
5. Treatment plan.....	15
6. What to expect with Silk'n Model H5003	16
7. Get started!.....	17
8. Cleaning & maintenance	19
9. Parts and consumables	20
10. Troubleshooting	21
11. Customer service	21
12. Warranty	22
13. Disposal	22
14. Specifications	23
15. Signs & symbols	25
16. Labeling	26
17. Compliance Statements	27

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Contraindications

- ⚠ Silk'n Model H5003 must not be used by children, adults with limited physical, sensory or psychological capacities, or those lacking experience or knowledge in how to use the device. Keep the device out of the reach of children.
- ⚠ Use a lower EMS energy level if you experience any discomfort.
- ⚠ Stop use immediately if your skin is very red, or if it blisters or burns.
- ⚠ Use on clean, dry, healthy skin.

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2. Safety

With Silk'n safety comes first. Silk'n Model H5003 has these safety features:

- **Proximity sensor:** The device works only when it is on, and the treatment unit/s are in good contact with your skin.
- **Auto shut-down:** Although you can stop a treatment whenever you wish, the device shuts down after a minute or so, when there is no contact with the skin (the *operating light* will be constant green). It also shut downs automatically after 15 minutes of continuous use.

3. Get to know Silk'n Model H5003

Silk'n Model H5003 is a system to help you exercise your muscles and tone specific parts of your body without physical effort, and in the comfort of your home. Target areas are the upper arms, abdomen, buttocks or thighs.

Silk'n Model H5003 is a powerful electrical device and should be used with special attention to safety. **Before use, it is important that you read and understand the instructions on how to use the device including after-use procedures, and follow them strictly.** We recommend you re-familiarize yourself with the instructions in this manual before each treatment.

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3.2 Device description

Silk'n Model H5003 is a device for home use that utilizes electrical muscle stimulation (EMS). This technology helps increase muscle tone and enhances body toning and reshaping.

Silk'n Model H5003 includes a selection of belts suited to the part of the body you wish to treat. Each treatment unit contains electrodes that activate the EMS. The treatment unit/s are attached to the belt and then placed on placed on the part of the body you wish to tone and reshape. Silk'n Model H5003 uses rechargeable batteries so the treatment unit/s can be easily charged before each use.

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3.3 Intended use

Silk'n Model H5003 is a non-invasive aesthetic treatment to improve muscle tone and firmness, and strengthen muscles in the waist, arms, thighs and buttocks.

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
4. Device controls & settings

4.1 Connecting the different parts


Connecting Silk'n Model H5003 is easy. First connect the *treatment unit/s* to the *belt*, and then wear the *belt* on the part of the body you wish to treat.

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

7. Of course, to use Silk'n Model H5003, you will need to disconnect one side, place the *belt* on your body and then re-connect the side you disconnected!

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4.2 Controls and settings

Battery / Error Operating Not in use

Each *treatment unit* has a control panel (#2 and #3 in diagram at the front of the manual). The control panel includes:

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4.3 Energy levels

Silk'n Model H5003 has 15 levels of Electrical Muscle Stimulation (EMS). This gives you excellent control of the intensity of muscle stimulation you wish to use for each part of the body you treat. What's more, you can change levels as often as you like during a treatment.

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

1. Connect the empty *stand/charging cradle* to the electricity supply. The 2 *charging indicator lights* will be constant green.
2. Place the *treatment unit/s* in the cradle so that they sit in the designated grooves.
 - When the *treatment unit* is charging, the orange *battery/error indicator light* on the *treatment unit* will flash on and off.
 - When fully charged the *battery/error indicator light* (left) will go off, and the green *operating light* (middle) will be constant green.

3. Once charged, disconnect the *stand/charging cradle* from the electricity supply.

5. Treatment plan

EN

- ⚠ For use on clean, healthy skin (see section "1. Warnings & contraindications").
- Treat the target area for 15 minutes per treatment session.
- For larger parts of the body, like the stomach or buttocks, you may wish to move the *treatment unit/s* to an adjacent part of the target area, and use for another 15 minutes. You may keep doing this until you have covered the target area.
- Try to avoid gaps and overlaps.
- Never use Silk'n Model H5003 on exactly the same spot during a treatment session.
- For daily use, but no more than 3 times a week on each specific target area. Continue until you get the results you want.

6. What to expect with Silk'n Model H5003

If you follow the instructions fully and pay attention to **all warnings and contraindications**, complications and side effects are uncommon. However every energy-based device, including those designed for home use, involves some degree of risk. Therefore it is important that you understand and accept the risks and complications that may occur. Be sure to review section "1. Warnings & contraindications".

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7. Get started!

You're on your way to a trim and toned body! Treating with Silk'n Model H5003 is easy. Just follow these steps for each part of the body you wish to treat.

7.1 Treatment steps

1. Use on clean, healthy skin (refer to section "1. Warnings & contraindications").
 2. Check that the batteries are well charged before use. If the *battery/error light* is orange and blinks slowly, this indicates that the battery is low.
 3. Before each treatment, wet the *electrode pad/s* with a wet cloth or with wet fingers, and place them on each *treatment unit* (see section "4.4 Wetting the electrode pads"). Next, connect one or more *treatment units* to the *belt* (see "4.1 Connecting the different parts").
 4. Disconnect one of the *connectors* on one side, and then place the *belt* around the area you wish to treat. Then reconnect the side you just disconnected.
 5. Adjust the tightness of the *belt* so that the *treatment unit/s* are in good contact with your skin, and that the *belt* is firm and well secured to your body.
-  If the *belt* is not tight against your skin, the EMS energy may disconnect. If that happens, press the + *button* to reselect the EMS energy level required.
6. **Switch on.** With the controls in view, press the *control button* on each *treatment unit*. The green and orange *indicator lights* will go on, and then go off after a second or so and the *treatment unit/s* will vibrate once. When the *treatment unit/s* are in good contact with the skin, the green *operating light* will blink slowly.

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External: Nylon. Middle
layer: SBR. Internal: Nylon

8.5 Storage

- After cleaning, we recommend that you return the *treatment unit/s* to the *stand/charging cradle*.
- If you do not plan to use the device for some time, store the different parts in the original packaging in a clean, dry place to prevent the accumulation of dust.
- Make sure the *stand/charging cradle* is disconnected from the electricity supply.

9. Parts and consumables

You can purchase replacement treatment unit/s or belts through the retailer where you bought the device, or you can order directly from www.silkn.com.

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14. Specifications

Model No.	H5003 (device); H5002 (charging cradle)
Technology	Electrical Muscle Stimulation (EMS)
Charging cradle	Input 12Vdc, 1.5Amax Charging 2 modules simultaneously
Module powered by	Li-ion battery, 3.7V, 2600mA
Adapter models	PMA18A-120 (USA) KSA24A1200150HU (USA)
Adapter Rating Input: Output:	100-240 Vac, 50-60 Hz, 0.6A Max. 12 Vdc, 1.5A 12 Vdc, 1.5A
EMS	0-32mA into 1000 ohm load
- Pulse Amplitude	
- Output Voltage	Adjustable 0-32V, Max output 32V peak to peak into 1000ohm load.
- On Time	Adjustable, 20~2400 milliseconds
- Off Time	Adjustable, 20~2400 milliseconds
- Ramp Time	Adjustable, 22~800 milliseconds, 22 milliseconds/step. The "On" time will increase and decrease.
- Pulse Rate	30~50 Hz
- Pulse Width	350~450µs
- Timer	15 min
- Wave Form	Symmetrical Bi-phasic rectangular pulse
- Tolerance	There may be a +/-5% tolerance of all settings and +/- 20% tolerance of output of intensity.
Package size	192 x 296 x130 (L x W x H) [mm]

Weight	Module: 300 [gr] Charging Cradle: 215 [gr]
Transport & storage between uses and storage condition	Temperature: minus 40 - 70 [°C] Relative humidity: 10% - 90% Atmospheric pressure: 500 - 1060 [hPa]
Operating conditions	Temperature: 5 - 40 [°C] Relative humidity: 15% - 90% Atmospheric pressure: 700 - 1060 [hPa]
Time from minimum storage temperature between uses to minimum operating temperature with ambient temperature of 20°C degrees	30 minutes
Time from maximum storage temperature between uses to maximum operating temperature with ambient temperature of 20°C degrees	30 minutes
Service life	3 years

15. Signs & symbols

EN

	Declaration of Conformity: Products labeled with this symbol meet all applicable provisions of the European Economic Area.
	CSA mark for US and Canada.
	Assembled from tested components.
IP22	The device is protected against the ingress of objects such as fingers.
	Degree of protection against electric shock. BF applied part.
	Keep dry.
	Do not dispose with household waste! The WEEE directive requires that you dispose of this electrical device (including all accessories and parts) at a designated recycling collection point for electrical devices.
	Follow the operating instructions.

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(01)07290016027600(11)171030(21)L11744000001

This device complies with Part 15 of the FCC Rules.
Operation is subject to the following two conditions:
(1) This device may not cause harmful interference
and (2) This device must accept any interference
received, including interference that may cause
undesired operation.

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Table 2: Guidance and Manufacturer's Declaration: Electromagnetic Immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the devices should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 15 kV air	± 8 kV contact ± 15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV for power supply lines ± 1 kV for input/output lines	± 2 kV for power supply lines NA	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	2 kV 0, 90, 180 & 270 degree phase angles	2 kV 0, 90, 180 & 270 degree phase angles	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	100% drop, 0.5 periods (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) 100% dip, 1 period 30% dip, 25/30 periods 100% drop, 5 seconds	100% drop, 0.5 periods (0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°) 100% dip, 1 period 30% dip, 25/30 periods 100% drop, 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note UT is the a.c. mains voltage prior to application of the test level.			

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Immunity Test	IEC 60601-1-2 Test Level	Compliance Level	Electromagnetic Environment Guidance
Conducted RF IEC 61000-4-6	6 Vrms 150 kHz to 80 MHz	6 Vrms 150 kHz to 80 MHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance. $d = \sqrt{\frac{3S}{P_1}}$
Radiated RF	10 V/m	10 V/m	$d = \sqrt{\frac{3S}{E_1}}$ 80 MHz to 800 MHz
IEC 61000-4-3	80 MHz to 2,7 GHz	80 MHz to 2,7 GHz	$d = \sqrt{\frac{7}{E_1} N P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 

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Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz $d = \left[\frac{3,5}{V_1}\right]\sqrt{P}$	80 MHz to 800 MHz $d = \left[\frac{3,5}{E_1}\right]\sqrt{P}$	800 MHz to 2,5 GHz $d = \left[\frac{7}{E_1}\right]\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.7	3.7	7.38
100	11.7	11.7	23.33
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. Note 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

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A distance of at least 0.5cm. between the equipment's antenna and all persons should be maintained during the operation of the equipment.

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