

What's Included

- iHelp Max
- AC Adapter
- Charging Cable
- Breakaway Lanyard

Optional Accessories

- Charging Base
- Snap-On Belt Clip
- Wristband

Specifications

Construction

- Comfortable to wear
- Lightweight - 2 ounces
- Small size - 2.8" x 2.1" x .8"
- Adjustable breakaway neck lanyard
- Water-resistant in the shower (IP67)
- Durable and shatter resistant

Care & Maintenance

To extend the life of your iHelp Max, try to:

- Avoid using and storing it in dusty places
- Avoid leaving it in overheated or overcooled places
- Clean it regularly with a dry cloth
- Don't use chemicals or detergents
- Avoid disassembling or refitting the device

Ways to Wear your iHelp Max

Your iHelp Max device should be worn at all times. You can wear it around your neck with the included lanyard, simply carry it in a pocket or use one of our optional accessories to wear it on your wrist or belt.

Attaching the Lanyard

Our lanyard is designed with a breakaway safety feature to help avoid accidents.

Please Note:

By using any type of lanyard, you run the risk of choking, including the possibility of serious injury. The lanyard can get entangled in wheelchairs, walkers, bed guard rails, jewelry, clothes, or other objects. Please exercise caution when using this accessory.