

## ALL IMPLIED WARRANTIES ARE LIMITED TO THE PERIOD OF THIS LIMITED WARRANTY.

This warranty gives you specific legal rights and you may also have other rights which vary from state to state. Some states do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation or exclusion may not apply to you.

Every Ovutek product automatically includes a 2-year warranty.  
This warranty is made by: Ovutek

## Warranty & Customer Support

Your product comes with a 2-year warranty starting from the date of purchase. If you have any questions or concerns about your new product, please contact our helpful Customer Support Team.  
In order to solve your problem more quickly, please have your order invoice and order ID ready before contacting Customer Support.

Email: [service@ovutek.com](mailto:service@ovutek.com)

Manufacturer: Ovutek

Distributed by Ovutek

## Online Manual



Need more help?  
Scan here to view the manual for more information.

You can also type the following link in the web browser:  
<https://www.ovutek.com/manual/index?model=VFS335>


# Ovutěk™

**Royyt+ Mini Body Fat Scale**  
Model: VFS335-PUS


## Get Started




1. Download the Ovutek FitPal app
2. Create an account
3. Go to "Device Management"  
And follow device addition Guidelines



Tap the "+" button in the  
Ovutek app



Ensure the Bluetooth of  
your phone is on



Press the scale to activate  
the screen. Select the  
device in the app

## Using Your Scale

1. Place the scale on a hard, flat surface ( not a carpet or mat ).
2. Just step on, the scale will show your weight in real time.
3. The display will flash several times before showing your final measurement.






## Tips

1. Before first use, open the battery cover and peel off the insulation film.
2. To switch the unit of measurement ( kg, lb ), press the UNIT button on the under-side of the scale or adjust the settings in the app.
3. To prevent damage, please avoid hitting the edges of the scale against hard objects or surface.
4. **SLIPPERY WHEN WET!** Do not step on scale with wet feet.
5. If you have any questions or concerns about your new product, please contact our helpful Customer Support Team.

## Specifications

Model	VFS335-PUS
Weight Capacity	400 lb / 180 kg
Increments	D = 0.2 lb / 0.1 kg
Weight Units	lb / kg
Platform	Tempered Glass
Power Supply	2 x 1.5V AAA Batteries
Auto-Off	15 Seconds
Dimensions	10.2 x 10.2 x 0.9 in / 26 x 26 x 2.3 cm
Weight	2.2 lb / 1.0 kg
Operating Systems	iOS® 12.0+ / Android™ 7.0+

## Display Icons

	Low Battery (Need change batteries)
	Connection enabled for phone app syncing
	Zero-Current Mode is activated ( Electric current disabled. Scale will only measure weight and BMI ).
	The result maybe inaccurate due to some reason, we recommend you re-measure.
	Overload. Weight exceeds 403 lb / 183 kg

## Warranty Information

Product Name:	Royyt+ Mini Body Fat Scale
Model:	
For your own reference, we strongly recommend that you record your order ID and date of purchase.	
Order ID:	
Date of Purchase:	

## Terms & Policy

Ovutek warrants this product to the original purchaser to be free from defects in material and workmanship, under normal use and conditions, for a period of 2 years from the date of original purchase.

Ovutek agrees, at our option during the warranty period, to repair any defect in material or workmanship or furnish an equal product in exchange without charge, subject to verification of the defect or malfunction and proof of the date of purchase.

**There is no other express warranty. This warranty does not apply:**

- If the product has been modified from its original condition;
- If the product has not been used in accordance with directions and instructions in the user manual;
- To damages or defects caused by accident, abuse, misuse or improper or inadequate maintenance;
- To damages or defects caused by service or repair of the product performed by an unauthorized service provider or by anyone other than Ovutek;
- To damages or defects occurring during commercial use, rental use, or any use for which the product is not intended;
- To damages or defects exceeding the cost of the product.

Ovutek will not be liable for indirect, incidental, or consequential damages in connection with the use of the product covered by this warranty.

This warranty extends only to the original consumer purchaser of the product and is not transferable to any subsequent owner of the product regardless of whether the product is transferred during the specified term of the warranty.

## **FCC Statement**

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.