




Product: Wireless Probe Type
Ultrasound Scanner

Model: 5CB FCC ID : 2A3TY-5CB

Built-in Battery: 3.8 V  / 4200 mAh



 **Manufactured by Bionet Co., Ltd.**

5F, 61 Digital-ro 31-gil Guro-gu,
Seoul 08375, REPUBLIC OF KOREA

TEL : +82-2-6292-6410

FAX : +82-2-6499-7789


Operation Mode of operation: 1:2 duty cycle
(ON: 5 min/OFF: 10 min)

This device complies with part 15 of the
FCC Rules. Operation is subject to the
following two conditions:


(1) this device may not cause harmful
interference, and (2) this device must
accept any interference received,
including interference that may cause
undesired operation.

Product: Wireless Probe Type
Ultrasound Scanner

Model: 10L FCC ID : 2A3TY-5CB

Built-in Battery: 3.8 V  / 4200 mAh



 **Manufactured by Bionet Co., Ltd.**

5F, 61 Digital-ro 31-gil Guro-gu,
Seoul 08375, REPUBLIC OF KOREA

TEL : +82-2-6292-6410

FAX : +82-2-6499-7789


Operation Mode of operation: 1:2 duty cycle
(ON: 5 min/OFF: 10 min)

This device complies with part 15 of the
FCC Rules. Operation is subject to the
following two conditions:


(1) this device may not cause harmful
interference, and (2) this device must
accept any interference received,
including interference that may cause
undesired operation.

Product: Wireless Probe Type
Ultrasound Scanner

Model: 10LB FCC ID : 2A3TY-5CB

Built-in Battery: 3.8 V  / 4200 mAh



 **Manufactured by Bionet Co., Ltd.**

5F, 61 Digital-ro 31-gil Guro-gu,
Seoul 08375, REPUBLIC OF KOREA

TEL : +82-2-6292-6410

FAX : +82-2-6499-7789

Operation

Mode of operation: 1:2 duty cycle
(ON: 5 min/OFF: 10 min)

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

label location



ID WAPCAMA022