

# **RF Exposure Evaluation Report**

Product Name : Automatic Upper Arm Type Blood Pressure Monitor

Model No. : HL858DC

FCC ID : 2ABTAHNL85DC1

Applicant : Health & Life Co., Ltd.

Address : 9F, No. 186, Jian Yi Road, Zhonghe District, New Taipei City, Taiwan

Date of Receipt :	Apr. 09, 2021
Date of Declaration :	Apr. 23, 2021
Report No. :	2140228R-E3082100014
Report Version :	V1.0



The test results relate only to the samples tested.

The test results shown in the test report are traceable to the national/international standard through the calibration report of the equipment and evaluated measurement uncertainty herein.

This report must not be used to claim product endorsement by TAF or any agency of the government.

The test report shall not be reproduced without the written approval of DEKRA Testing and Certification Co., Ltd. Measurement uncertainties evaluated for each testing system and associated connections are given here to provide the system information for reference. Compliance determinations do not take into account measurement uncertainties for each testing system, but are based on the results of the compliance measurement.



Issued Date: Apr. 23, 2021 Report No.: 2140228R-E3082100014

# DEKRA

Automatic Upper Arm Type Blood Pressure Monitor				
Health & Life Co., Ltd.				
9F, No. 186, Jian Yi Road, Zhonghe District, New Taipei City, Taiwan				
Health & Life Co., Ltd.				
Living Science Co., Ltd.				
Health & Life Medical Technology Company Limited				
HL858DC				
2ABTAHNL85DC1				
Health & Life				
KDB 447498 D01 v06	<ul> <li>☐ Minimum test separation distance ≥ 20 cm</li> <li>☑ For low power devices</li> </ul>			
Complied				
:	peggy Tu			
	(Adm. Assistant / Peggy Tu )			
:	Wenlee (Senior Engineer (Wenlee)			
( Senior Engineer / Wen Lee ) : (Director / Vincent Lin )				
	Health & Life Co., Ltd. 9F, No. 186, Jian Yi Road Health & Life Co., Ltd. Living Science Co., Ltd. Health & Life Medical To HL858DC 2ABTAHNL85DC1 Health & Life KDB 447498 D01 v06			



# 1. GENERAL INFORMATION

## **1.1. EUT Description**

Product Name	Automatic Upper Arm Type Blood Pressure Monitor	
Trade Name	Health & Life	
Model No.	HL858DC	
FCC ID.	2ABTAHNL85DC1	
Frequency Range	2402-2480MHz	
Channel Number	V4.2: 40CH	
Type of Modulation	V4.2: GFSK(1Mbps)	
Antenna Type	PCB Antenna	
Channel Control	Auto	
Antenna Gain	Refer to the table "Antenna List"	

#### Antenna List

No.	Manufacturer	Part No.	Antenna Type	Peak Gain
1	atl	LB-BLE-005	PCB Antenna	2.23dBi for 2.4 GHz

### 2. **RF Exposure Evaluation**

#### 2.1. Standard Applicable

According to 1.1307 (b)(1), systems operating under the provisions of this section shall be operated in a manner that ensure that the public is not exposed to radio frequency energy level in excess of the Commission's guideline.

#### 2.2. Measurement Result:

According to KDB Publication 447498 D01, section 4.3.1, per the calculations of item 1 (Power(mW)/separation (mm)\*sqrt(f(GHz) $\leq$ 3.0), SAR is required as shown in the table below where calculated values are greater than 3.0:

1.) Operation frequency = 2450MHz and antenna separation distance = 5mm, SAR Test Exclusion Threshold = 10mW

	Maximum output power			SAR Test	
Frequency Band	Peak Gain: 2.23dBi		Exclusion Threshold	Calculated Threshold Value	
(MHz)	Power	EIRP	EIRP	(mW)	( $\leq$ 3.0 SAR is not required)
	(dBm)	(dBm)	(mW)	(mW)	
2402 - 2480	4.29	6.52	4.49	10	1.402

Note1: The SAR/MPE measurement is not necessary.

Note2: The maximum peak output power is refer to report No.: 2140228R-E3032110108 from the DEKRA.