

1. PLACE ITEM IN POUCH.
2. PEEL OFF PROTECTIVE STRIP FROM ADHESIVE BAND.
3. FOLD OVER PAPER FLAP ALONG PERFORATION AND PRESS DOWN FIRMLY FROM THE CENTRE OUTWARDS.

Sterilization Date:

/ /

Stellar Telemetry

Model # 43001-IMP-02

FCC ID: 2AC4C-AU430001SM

IC: 12302A-AU430001SM



Made in USA



Quantity: P23-0149

3 1/4" x 6 1/2" (83mm x 159mm)

Processed yellow

Don't use if package is damaged. ISO11607-1
Processed dark grey