

Mark Bertinetti  
24 July 2014

**Common labels/moulded in artwork**




**BRAND**  
ResMed branding  
chrome or colour foiled  
on front of product




**“BF APPLIED PART”  
Body Fibrillation Icon  
& Text**  
Moulded into plastic case  
using an insert  
near rear Outlet



**“REFER TO MANUAL” Label  
Man reading book Icon & Text**  
White pre-printed label stuck next to  
power connector on rear of chassis



**REGULATORY Label**  
White pre-printed label  
stuck on bottom of chassis




**FCC & IC Compliance Labels**  
White pre-printed label  
stuck on bottom of chassis

**Product Specific Labels**




**37015 AirSense10 CPAP - USA**

R370-7204



**37020 AirSense10 Elite - USA**

R370-7205






**37028 AirSense10 AutoSet - USA**

R370-7206



**37036 AirSense10 AutoSet FH - USA**

R370-7207



**37043 AirSense10 ASV - USA**

R370-7210






**37046 AirSense10 S - USA**

R370-7208



**37051 AirSense10 VAuto - USA**

R370-7209



**PRODUCT NAME PLAQUES**  
Foiled product variant badges  
“permanently” attached to front of product

**SERIAL NUMBER LABELS**  
Thermally printed labels  
attached to rear of product



**Figure 1 Front of CPAP**



Figure 2 Back of CPAP



**Figure 3 Underneath of CPAP**