

ResMed Confidential

Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way. © 2022 ResMed.

SPEC NO.	eLB1392-074
REV	6.0
PAGE	1 of 4

Melanie Theaker 8 September 2022

Common labels/moulded in artwork



BRAND


ResMed branding chrome or colour foiled on front of product



"BF APPLIED PART"
Body Fibrillation Icon
Moulded into plastic case
using an insert
on rear of device,
adjacent to the device outlet




"R390-411 Label - MR Unsafe"
White pre-printed label stuck next to power connector on rear of assembly



"R390-429 Label - FCC & IC NIR M1"

White pre-printed label stuck on bottom of assembly





A390-7374 Bottom Case Assembly Laser-Etched on bottom of chassis



ResMed Confidential

Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way. © 2022 ResMed.

SPEC NO.	eLB1392-074
REV	6.0
PAGE	2 OF 4

Melanie Theaker 8 September 2022

Product Specific Labels for FCC







NOTE:

- The device ID is the base code, which can be packaged in different boxes of accessories
- CO means it is packed with a humidifier.
- TRI means it is packed with a humidifier and heated tube.



ResMed Confidential

Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way. © 2022 ResMed.

SPEC NO.	eLB1392-074
REV	6.0
PAGE	3 of 4

Melanie Theaker 8 September 2022

Product Specific Labels for ISED







NOTE:

- The device ID is the base code, which can be packaged in different boxes of accessories
- CO means it is packed with a humidifier.
- TRI means it is packed with a humidifier and heated tube.



ResMed Confidential

Proprietary Information, not to be reproduced or made available to third parties without prior consent from ResMed and not to be used in any unauthorised way. © 2022 ResMed.

SPEC NO.	eLB1392-074
REV	6.0
PAGE	4 OF 4


Melanie Theaker 8 September 2022



REAR VIEW OF DEVICE



BOTTOM VIEW OF DEVICE



Additional Label for AMER devices fits into cavity