	<i>Labelling</i>	<i>Reference</i>
	Plan Labelling ePerf Connect laser label	P110LAB003

	Name	Position	Date	Signature
Writing	Hélène Bataille	Assistante marketing et commerciale		
Verification	Liang ZHOU	Engineer R&D		
Approval	BLOND Estelle	Technical Director		

Modification History

Revision	Modification Date	Modifications
1	11/03/2021	Creation
2	12/08/2021	Add BF symbol, modification the size of QR code and SN symbol and fond
3	18/01/2022	Modification QR code size and add manufacturing date in the QR code Add USB, Model, input information
4	01/12/2022	ADD FCCCID and move the BF symbol's position
5	09/12/2022	ADD FCCID for BLE module and change the text's size

Page intentionally left blank

1. Purpose of the document and scope of application

This document describes the laser label plan for ePerf Connect.

2. Documents and acronyms

A. Applicable documents

NA

B. Reference documents

PRO18: Plan, Design and Develop the product.

C. List of acronyms and definitions

NA

Label Information :

Support : Laser
 SN: Serial Number

The QRcode includes:

BodyCAP GATEWAY

Ref : P110


SN : XXXX

Manufacturing date : XXXX

<https://www.bodycap-medical.com/>

Position on the product :

This label is printed on the bottom shell of the e-Celsius medical connect system



Real size
 1:1

