


## Label & Location

Back Side



**MEDICAL-GENERAL MEDICAL EQUIPMENT AS TO ELECTRICAL SHOCK,  
FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH  
ANSI/AAMI ES 60601-1(2005)+AMD(2012) and CAN/CSA-C22.2**

**No. 60601-1 (2008) + (2014)**

**Product Name: Tablet PC**

**Model : CyberMed Rx**

**Power : 19Vdc  $\approx$  3.78A**

**FCC ID: 2AHZW-CYBRXX10  
IC: 24472-CYBRXX10**



**4WW4  
E357129**

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

**USE EDAC POWER ELECTRONICS CO., LTD**

**Adapter model EM10681V, AC Input 100-240 Vac, 2.0-1.0 A, 50-60 Hz**

**Tested To comply With FCC class B Standards**

**www.cybernet.us  
RoHS/WEEE**



**Product Name: Tablet PC**

**Power : 19Vdc  $\approx$  3.42A**

**Model : Rugged X10**

**USE FSP GROUP INC.**

**Adapter Model : FPS065-REBN2 , AC Input 100-240 Vac, 1.5A, 50/60Hz**

**Tested To comply With FCC class B Standards**

**This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.**

**RoHS/WEEE**



**www.cybernet.us**

**FCC ID: 2AHZW-CYBRXX10  
IC: 24472-CYBRXX10**

