


Theraflu Flu Test (I-FLU-C02)



Details of decal label (sticker): position, content and dimensions

Position on EUT:


Blue 'false coloring' indicate position of the decal on back cover of device



Label information (content):

Font: Arial


Size: 4.5



Label application method: manual

Label material: Polypropylene (Sticker)

Decal dimensions



Wireless Analyser Compliance Decal
Scale 1:1

Do not print background wash,
dimensions or knifelines for visual guide only.



GENERAL NOTES

- 100 Micron
- White Vinyl
- Adhesive Backed



Formzoo Design Pty Ltd
PO Box 110A, Spring Hill, Queensland, Australia 4008
Mobile: +61 7 3268 7228 Email: info@formzoo.com.au
Tel: +61 7 3461 5475 URL: www.formzoo.com.au

client: Ellume Pty Ltd

project: GSK Home Flu Test

job No: zoo615

Copyright 2019
This drawing is supplied in confidence and should not be used for any
purpose other than that specified. Do not disseminate this drawing
without prior written consent from Formzoo Design Pty Ltd.

issue: 02 DATE: 25-06-2019

DESCRIPTION: Wireless Analyser Compliance Decal

PARTY: zoo615-g507