

E-labeling on HUB01 screen


The “FCC ID: 2A08L-HUB01 contains FCC ID: 2A08L-WL18DBMOD” text is included in the Portrait Mobile Patient Monitor HUB01 e-labeling as presented in the figure below. The e-labeling is visible for an end user via the following steps:

- Press the physical HUB01 “star button”
- → Choose “Regulatory Markings” from the HUB01 touch screen GUI
- → “Regulatory Markings” page (below) opens to the HUB01 screen

Portrait HUB01, Mobile Patient Monitor


HUB01 GUI touch screen

HUB01 UI star button




Physical labeling in HUB01 device

The “FCC ID: 2A08L-HUB01” text is laser etched to the HUB01 back cover as presented in the figure below. The font used is “GE Inspira Sans”, the font size is 6pt and the font color is “Pantone: Cool Gray 10”



Physical labeling in HUB01 sales package

The "FCC ID: 2A08L-HUB01" text is printed to the HUB01 sales package as presented in the figure below.



XXXXXXXXXX replaced with the text 2A08L-HUB01