

Begin Your Session

- Customize the temperature and dose to your desired preferences.
- Keep track of
 - How many draws you've had
 - How much you've inhaled per draw
 - The total amount inhaled
 - The average temperature of your session



Record How You Feel

- Rate how much you enjoyed the session
- Select any mental or physical side effects you experienced during or after your session
- Track how effective your session was at relieving the symptoms of your medical condition



Cleaning Mouthpiece

- Step 1:**
- Detach the mouthpiece by pressing upwards on the front edge




Cleaning Heating Chamber

- Step 1:**
- Before cleaning, make sure heating chamber temperature indicator is GREEN before proceeding
 - If temperature indicator is OFF, please turn on RYAH and ensure temperature indicator is GREEN before proceeding



Caution: RED temperature indicator means chamber is hot!

- Step 2:**
- Pull out the bottom chamber lid



- Step 3:**
- Soak the mouthpiece and extended pathway / insert in soapy water or isopropyl alcohol for 15 minutes. Do not use acetone or bleach as they can damage or more than 15 minutes could affect exterior finish!

- Step 4:**
- After soaking for 15 minutes, use the scraping tool to remove any extra residue from the mouthpiece and extended vapor pathway

- Step 3:**
- Detach mouthpiece by pressing upwards on the front edge



- Step 4:**
- Carefully insert the pipe cleaners included in the RYAH cleaning kit through the bottom chamber and mouthpiece openings to remove any residue within the heating chamber.



WHAT'S IN THE BOX
 RYAH Device
 Charger
 Cleaning Kit
 Manual