























Model name: DVM-A1

		Vathin® Digital Video Monitor	
ES:Monitor de video digital		DE:Digitaler Videomonitor	
FR:Moniteur vidéo numérique		IT:Monitor video digitale	
PT:Monitor de video digital		KR: 디지털 비디오 모니터	
DVM-A1			
REF XXXXXXXXX			
Input:15 VDC, 4.2A max			
Internal: 10.8 VDC, 9000mAh			
 XXXX-XX-XX	SN QXXXXXXX	Rx Only	
FCC ID:2AY4E-DVM	IC :27001-DVM		
5150-5250MHz is only for indoor use.			
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.			
 Hunan Vathin Medical Instrument Co., Ltd. 1/F, Building 12, Innovation and Entrepreneurship Service Center, No. 9 Chuanqi west road, Jiuhua Economic Development Zone, Xiangtan, Hunan, China			
Lotus NL B.V.			
Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands. peter@lotusnl.com +31645171879, +31626669008		(01)0 6972197480566 (11)200318 (21)Q2003002	
EC	REP		
			
(01) 0 6972197 48056 6 (11) 200318 (21) Q2003002			

Model name: DVM-A2

		Vathin® Digital Video Monitor	
ES:Monitor de video digital		DE:Digitaler Videomonitor	
FR:Moniteur vidéo numérique		IT:Monitor video digitale	
PT:Monitor de video digital		KR: 디지털 비디오 모니터	
DVM-A2			
REF XXXXXXXXX			
Input:15 VDC, 4.2A max			
Internal: 10.8 VDC, 9000mAh			
 XXXX-XX-XX	SN QXXXXXXX	Rx Only	
FCC ID:2AY4E-DVM	IC :27001-DVM		
5150-5250MHz is only for indoor use.			
This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.			
 Hunan Vathin Medical Instrument Co., Ltd. 1/F, Building 12, Innovation and Entrepreneurship Service Center, No. 9 Chuanqi west road, Jiuhua Economic Development Zone, Xiangtan, Hunan, China			
Lotus NL B.V.			
Koningin Julianaplein 10, 1e Verd, 2595AA, The Hague, Netherlands. peter@lotusnl.com +31645171879, +31626669008		(01)0 6972197480573 (11)200318 (21)Q2003002	
EC	REP		
			
(01) 0 6972197 48057 3 (11) 200318 (21) Q2003002			

Label Location

