

Label & Location



Label position

Label

AMSINO

AMSure® Enteral Feeding Pump **REF**ASEP00

SN XXXXXXXX

Model: ASEP00

 XXXX-XX

Adapter Input:  100-240V, 50/60Hz

Adapter Output:  15V

Adapter Model: LXCP26-150160A

Internal Power Supply: DC 11.1V



FCC ID: 2AYKB-ASEP00

This device complies with Part 15 of the FCC Rules.


Operation is subject to the following two conditions:

- (1) this device may not cause harmful interference, and
- (2) this device must accept any interference received,

Including interference that may cause undesired operation.



IPX4



Caution



Rx
ONLY

Please read User Manual carefully before use.

Please use AMSINO recommended or approved feeding sets only.

Operation steps: Install feeding pump→Power on→Install feeding set→RFID Identification→Parameters setting→Priming→Connect to gastric tube→Start feeding→Feeding complete→Power off.



Amsino International, Inc.

708 Corporate Center Drive

Pomona, CA 91768, USA

www.amsino.com/symbol-glossary

Made in China

L10076 Rev. A

