

Labels are placed on the bottom of the device, left and right hand side of a groove when looking at the bottom of the device.





Swatches

THERMAL TRANSFER PRINT

 PRINT - BLACK - VARIABLE / SERIALISED INFO
 PRINTER QUALITY - MINIMUM 600DPI


Key

-  Part outline
Reference only. Do not print.
-  Dimensions
Reference only. Do not print.

Notes


1. Material - Flexcon THERMLfilm SELECT 22720
2. Finish - N/A
3. Tolerances unless otherwise stated +/- 0.2mm
4. Variable data to be set at time of printing according to factory order
5. Static artwork proof, refer to page 2.
6. Supporting print files:
 PDF 1:1 Artwork
 Reference image 1500DPI transparent PNG
7. All variable data text to be GOTHAM - BOLD 6.5pt unless otherwise stated
8. Refer to page 2 for static artwork proof
9. Print Quality: To be printed on minimum 600DPI printer
10. Export: PDF 1:1 artwork, Reference image 1500dpi transparent PNG

REV	DATE	CO	DESCRIPTION	DRWN	APR
A	09-DEC-2020		RELEASE FOR EUA	TRG	BJD
B	17-DEC-2020		UPDATED LABEL ID	TRG	BJD
C	15-JAN-2021	CO 56229	BARCODE ADDED	TRG	BJD
D	27-JAN-2021	CO 56229	BARCODE POSITION CHANGED	TRG	BJD
E	25-MAR-2021	TBC	UPDATED REGULATORY INFO	TRG	BJD



SCALE 1:1

*8 x 8mm pink box
 Represents the 2d barcode position

 <p>Planet Innovation +61 3 9945 7510 planetinnovation.com.au</p>	© This drawing is confidential and must not be copied or disclosed in part or full without written consent of Planet Innovation	PI NO. PI-0074044	REV: E	
	TITLE: LEFT LABEL, BASE COVER, RADIATA READER EUA		SCALE: DO NOT SCALE	
	DRN: TRG APR: BJD	26-MAR-2021 26-MAR-2021	CHANGE ORDER TBC	SIZE SHEET A3 1 OF 2

PINE TREES HEALTH

REF

SN

IVD For In Vitro Diagnostic Use

**Emergency Use
Authorization only**

Rx only

FCC ID: 2AYYC-2232021

IC: 27017-2232021

HVIN: PTHR10





ART-L1-EUA-E

Swatches

THERMAL TRANSFER PRINT

 PRINT - BLACK - VARIABLE / SERIALISED INFO
 PRINTER QUALITY - MINIMUM 600DPI


Key

-  Part outline
Reference only. Do not print.
-  Dimensions
Reference only. Do not print.


Notes

1. Material - Flexcon THERMLfilm SELECT 22720
2. Finish - N/A
3. Tolerances unless otherwise stated +/- 0.2mm
4. Variable data to be set at time of printing according to factory order
5. Static artwork proof, refer to page 2.
6. Supporting print files:
PDF 1:1 Artwork
Reference image 1500DPI transparent PNG;
7. All variable data text to be GOTHAM - BOLD 6.5pt unless otherwise stated
8. Refer to page 2 for static artwork proof
9. Print Quality: To be printed on minimum 600DPI printer
10. Export: PDF 1:1 artwork, Reference image 1500dpi transparent PNG

REV	DATE	CO	DESCRIPTION	DRWN	APPR
A DRAFT 01	27-NOV-2020	TBC	INITIAL RELEASE FOR REVIEW	TRG	BJD
A	03-DEC-2020	TBC	RELEASE FOR INVESTIGATIONAL USE ONLY	TRG	BJD
B	03-DEC-2020	TBC	UPDATED RELEASE REMOVAL OF UL LOGO	TRG	BJD
C	17-DEC-2020	TBC	UPDATED LABEL ID FOR EUA RELEASE	TRG	BJD
D	26-MAR-2021	TBC	UPDATED REGULATORY INFO FOR EUA	TRG	BJD




SCALE 1:1

 Planet Innovation <small>+61 3 9945 7510 planetinnovation.com.au</small>	<small>© This drawing is confidential and must not be copied or disclosed in part or full without written consent of Planet Innovation</small>	<small>PI NO.</small> PI-0071704	<small>REV:</small> D	
	<small>TITLE:</small> RIGHT LABEL, BASE COVER, RADIATA READER		<small>SCALE:</small> DO NOT SCALE	
	<small>DRN:</small> TRG	26-MAR-2021	<small>CHANGE ORDER</small> TBC	<small>SIZE</small> A3
	<small>APR:</small> BJD	26-MAR-2021		<small>SHEET</small> 1 of 2

PTHealth Reader



Manufactured for:
Pine Trees Health
700 Main Street,
Cambridge,
MA 02139, USA



ART-L2-EUA-D