


PRODUCT MANUAL



Hongkong Tian Jie International Co., LTD

Scope of application

【1】

relieve eye fatigue and prevent myopia

【2】

beauty remove black, activate eye skin

【3】

improve sleep, inhibit neurasthenia

【4】

massage the temples to relieve headache

【5】

pulse vibration, dredge the eye meridian

【6】

far infrared hot compress, promote eye blood circulation

【7】

to restore ciliary muscle elasticity, improve eyesight

【8】

the built-in magnet promotes cell metabolism and activates cells

Intended for



Office workers sedentary, blood circulation is slow, the brain blood supply blocked, visual fatigue serious, long time toward the screen, easily lead to eye interference, visual acuity decline.



Stay up late the family often stay up late, the eye itself overload work, long time with the eye, dark circles are serious, easily lead to dry eyes, swelling pain, pain, visual fatigue increased, resulting in decreased vision.



Drivers long time, excessive use of eyes, stay up late will lead to excessive visual fatigue of the eye.



Bow heads, time to use the phone, and toward the phone screen, will lead to visual fatigue and decline.



Students long-term learning and use of eyes, will cause visual fatigue, vision loss.



The beauty of this product for the beauty of black eye circles have a very good effect.



Insomnia group prone to insomnia, sleep quality is not good.



A group of myopia or high myopia.




Maternity or postpartum group.



Pay no attention to eye care groups.

Use function description



1: Switch key

Switch the machine for 2 seconds long. Switch mode short by 1 second.

2: Bluetooth key



Long press Bluetooth to turn on / off and prompt tone.

3: Next key


Press one second to play the next piece of music. Long press 2 seconds, music pause.

4: Volume key

The volume size is divided into 4 files, and the boot defaults 2 files. You can press small to large 1 - 4 files. Each press 1 times, 1 files.



Five modes of intelligent function



Warm prompt: each mode defaults to 15 minutes automatic shutdown. A warning tone occurs when the battery voltage is too low.

Notice of use

- 1) : Please charge the product for 2 hours before use. The product won't be full until it leaves the factory.
- 2) : Prohibit charging side use, charging time should not be too long, too long charge will affect battery life.
- 3) : Please remove your glasses or contact lenses before use.
- 4) : Before use, please adjust the headband to the most comfortable position. The mask is aligned with the product.
- 5) : If the whole family is used together, it is recommended to use a clean soft paper pad every time.
- 6) : The use of the product by sitting or prone position, make the product more fit for the eye, the effect is better.
- 7) : After using this product, please close your eyes for 3 - 5 minutes, make eye acupuncture meridians to relax stretch, clearer vision.
- 8) : The eyes have just been operated on and are prohibited from being used. The doctor is under the guidance of the doctor.
- 9) : Before using this product, you need to connect the Bluetooth device. Please connect the Bluetooth device and use it. It is forbidden to walk or drive in the process of using.
- 10) : This product banned washing or explosion so as not to affect the electronic component Sai, service life.

After sales service and support

YAOEOAY intelligent eye massager product warranty regulations.

The after-sale service of this product is strictly based on the <<law of the People's Republic of China on the protection of consumers' rights and interests >>and the <<quality of domestic products of the people's Republic of China>>, and the after-sale Three Guarantees service is as follows:

1: Since your receipt within 7 days from the date of this product, "fault products listed in table" performance fault situation, determined via Hongkong Tian Jie International Limited Customer Service Technical Service Center, you can enjoy a free return or replacement service.

2: Since you sign 8 within 30 days from the date of this product, "fault products listed in table" performance fault situation, determined via Hongkong Tian Jie International Limited Customer Service Technical Service Center, you can enjoy a free replacement or repair service.

3: Since you sign within 12 months from the date of this product, "fault products listed in table" performance fault situation, determined via Hongkong Tian Jie International Limited Customer Service Technical Service Center to check, can enjoy free maintenance service.

Product performance fault table

Product name	Performance failure	Remarks
Eye massager	1: The product cannot boot	
	2: The product cannot be connected to Bluetooth or unstable	
	3: The product indicator is not on	
	4: The product cannot connect to the computer	Extended memory version before you can download
	5: The product cannot play music	
	6: Product massage air sac or abnormal air pressure	
	7: Abnormal motor vibration of the product	
	8: The solenoid valve does not work	
	9: Battery damaged or leaking	
	10: The heating sheet is broken or not hot	
	11: Damaged or damaged outer casing	Not man-made damage
	12: Headband abscission	Not man-made damage
	13: Cannot charge	
	14: Main board function program disorder	
	15: Button drop	
	16: Airbag breakage	
	17: Inadequate accessories	
	18: Other hardware problems cause functionality to be unavailable	

Non warranty regulation

- 1: Unauthorized maintenance, misuse, collision, omission, misuse, entry, accident, alteration, incorrect use of the product, or tear up, alter, mark, anti fake mark.
- 2: More than three packages warranty period has been exceeded.
- 3: Damage due to force majeure.
- 4: Failure to comply with performance faults listed in the product performance table.
- 5: Human factors lead to product failure in the product performance fault table.

Warranty card

Dear customer, this warranty card is an important document to apply for warranty service. Please keep it in good order with the seller. Thank you!

User information	User name		E-mail	
	Postal address			
	Contact number		Postcode	
Product information	Product name		Product barcode	
	Product model			
Agent information	Name			
	address			
	Contact number		Postcode	
	Date of sale		Invoice number	
Remarks				
This form is confirmed by seller's seal, and the warranty card is equivalent to "three packages of documents" for the products specified in the three pack".				

Item name: Intelligent eye protector

Item number: YEY-2.0

Battery capacity: 1100mAh

Rated voltage: 5V/1A

Accessories: complete machine, instruction manual, warranty card, certificate, flannelette bag.

Customer hotline: 4008358005

Changan town of Guangdong city of Dongguan province chongtou village Hongye Road No. 4.

Hongkong Tian Jie International Co., LTD

Dongguan Tian Jie Industrial Co., Ltd.

All technical information shall be subject to substantial performance. Hongkong Tian Jie International Co., Ltd. reserves the right of final interpretation.



客服热线：4008358005

东莞市长安镇涌头社区宏业中路4号

FCC Warning Statement

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment. This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- - Reorient or relocate the receiving antenna.
- - Increase the separation between the equipment and receiver.
- - Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- - Consult the dealer or an experienced radio/TV technician for help.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

The device has been evaluated to meet general RF exposure requirement.