

August 13, 2013

BABT FCB
Forsyth House,
Churchfield Road,
Walton-on-Thames,
Surrey,
KT12 2TD

Federal Communications Commission
445 12th Street, SW
Washington, DC 20554

Subject: 7000MHB Attestation Letter Re: Manual Update Requested by FCC

Dear Sirs or Madams,

We, Philips Healthcare, Home Monitoring, Lifeline Systems Inc., 111 Lawrence Street, Framingham, MA 01702, (508) 988-3002, hereby attest that our 7000MHB device manual will be updated to add the following (highlighted yellow) as requested by the FCC:

For best results, hold the Button up and away from your chest when you are speaking to Lifeline. Maintain a separation distance of 10 mm/ 0.394 in when operating the device in the held to face configuration, to meet the FCC RF Exposure safety requirements.

Sincerely,



Kimberly Cook

Senior Engineering Manager
Philips Healthcare
Home Monitoring
508-988-3002
kimberly.cook@philips.com