



Physio-Control, Inc.

11811 Willows Road NE, P.O. Box 97006, Redmond, WA 98052

Tel 425.867.4000 Toll-free 800.442.1142

[www.physio-control.com](http://www.physio-control.com)

May 23, 2012

Re: Antenna Information

To whom it may concern,

Per FCC § 15.203, the 'antenna' for the TrueCPR device are inductive loop antennas or coils located within the Backpad and Chestpad assemblies of the device. The Backpad coil is used for transmission purposes only and the Chestpad coils are used for receiving the transmission. No other external antenna can be used with the device. For reference, see the attached Discovery System Block Diagram.

Sincerely,


A handwritten signature in black ink, appearing to read "Matt Hedlund", is written over a light blue horizontal line.

Matt Hedlund  
Compliance Engineering  
Physio-Control, Inc.  
11811 Willows Road  
Redmond, WA 98073-9706



# Discovery System Block Diagram

Back Pad



Chest Pad

