

**Date:** 9/12/2013

**Re:** Model 6299 Latitude™ Rock Springs Point of Care (Consult)

**Communicator:** FCC ID Labels and Locations

The Guidant Corporation (a wholly owned subsidiary of Boston Scientific Corporation doing business as Boston Scientific Cardiac Rhythm Management, hereinafter referred to as BSC Model 6299 LATITUDE™ Rock Springs Point of Care (a.k.a., Consult) is an in-clinic or hospital monitor for use with inductive telemetry enabled BSC implantable pulse generators in the Cognis, Teligen, Progeny, Ingenio, Insignia, and Altrua devices.

The Model 6299 labels are shown in the following figures.

**Figure 1. BSC Model 6299 – Bottom Label (for reference only)**

**LATITUDE Consult™ Communicator**


For technical support, call: 1-800-CARDIAC  
(227-3422)






Federal (USA) law restricts this device to distribution and use by or on the lawful order of a physician.



**MODEL 6299**

Use only the power supply provided with this device:  
15 VDC, 1.2 A, AC adapter  
Globtek Model: GTM41061-1818-3.0

FCC ID: ESCCRM629913  
US: ESCMM00B6299  
Contains FCC ID: POOWML-C40

 **Manufacturer**  
Boston Scientific  
4100 Hamline Avenue North  
St. Paul, MN 55112-5798 USA

  
MR  
Unsafe
   
  
Non-ionizing  
electromagnetic  
radiation
   
  
Follow  
instructions  
for use
   
  
Do not  
throw in  
trash
   
  
Made in USA
   

  
P/N: XXXXXX-XXX  
  
**S/N: XXXXXX**  
15597-001

**Figure 2. BSC Model 6299 – Bottom Label Location (for reference only)**



**Figure 3. BSC Model 6299 – Rear Label and Location (for reference only)**



**Figure 4. BSC Model 6299 – Wand Bottom Label and Location (for reference only)**



Note: The wand bottom label is from a prototype device and is provided for reference only.