

DYNAGEN™ MINI ICD DR
Model D022



RA: IS-1
RV: DF4-LLHH

SN

Serial number

Use by



DYNAGEN™ MINI ICD DR Model D022



IS-1 BI DF4-LLHH RA RV	Uncoated device Do not use if package is damaged	Consult instructions for use STERILE EO Sterilized using ethylene oxide	Do not reuse Do not resterilize	Contents = Pulse generator + Torque wrench + Literature
---------------------------------	---	--	------------------------------------	---

Date of manufacture	Serial number	50°C Temperature limitation
Use by	Lot number	Dangerous voltage

Manufacturer
Boston Scientific
4100 Hamline Avenue North
St. Paul, MN 55112-5798 USA

This device may not interfere with stations operating in the 400.150 – 406.000 MHz band in the Meteorological Aids, Meteorological Satellite, and Earth Exploration Satellite Services and must accept any interference received, including interference that may cause undesired operation.

FCC ID: ESCCRMG17912
402 - 405 MHz
Made in Ireland

15474-001
EN US



Use by

SN

Serial number

DYNAGEN™ MINI ICD DR
Model D022



RA: IS-1
RV: DF4-LLHH

DYNAGEN™ MINI ICD DR Model D022



Implantable Cardiac Defibrillator DR

Characteristics as shipped

Tachy Mode Storage
Tachy Therapy Available Shock, ATP
Pacing Mode Storage
Pacing Therapy DDDR
Sensor Accelerometer
Pace/Sense Configuration RA: BI/BI RV: BI/BI

Symbols that may be used on packaging

Manufacturer	Serial number	Uncoated device
Date of manufacture	Reference number	Dangerous voltage
Consult instructions for use	Temperature limitation	Open here
Do not resterilize	Lot number	CRT-D RA, RV, LV
Do not use if package is damaged	Investigational use only	ICD RA, RV
Do not reuse	Place telemetry wand here	ICD RV
Use by	Sterilized using ethylene oxide	
Authorized representative in the European Community		

Refer to the System/Reference Guide for all programmable parameters and values.

If the container seal is broken, Boston Scientific will not accept contents for credit or restock.

Caution: Do not store near magnets or magnetic devices.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician trained or experienced in device implant and follow-up procedures.




DYNAGEN™ EL ICD DR
Model D152









RA: IS-1
RV: DF4-LLHH







Use by 


Serial number 

DYNAGEN™ EL ICD DR Model D152



IS-1 BI DF4-LLHH 	 Uncoated device  Consult instructions for use  Do not reuse  Do not use if package is damaged  Sterilized using ethylene oxide  Do not resterilize	 Contents = Pulse generator + Torque wrench + Literature
---	--	---


 Date of manufacture  Use by	 Serial number  Lot number	 Dangerous voltage  Temperature limitation
---	---	---

 Boston Scientific
4100 Hamline Avenue North
St. Paul, MN 55112-5798 USA
Manufacturer

This device may not interfere with stations operating in the 400.150 – 406.000 MHz band in the Meteorological Aids, Meteorological Satellite, and Earth Exploration Satellite Services and must accept any interference received, including interference that may cause undesired operation.

FCC ID: ESCCRMG17912
402 - 405 MHz
Made in Ireland

15470-001
EN US

Serial number 

Use by 

DYNAGEN™ EL ICD DR
Model D152

RA: IS-1
RV: DF4-LLHH


DYNAGEN™ EL ICD DR Model D152

























Implantable Cardiac Defibrillator DR

Characteristics as shipped

Symbols that may be used on packaging

Tachy Mode Storage
Tachy Therapy Available Shock, ATP
Pacing Mode Storage
Pacing Therapy DDDR
Sensor Accelerometer
Pace/Sense Configuration RA: BI/BI RV: BI/BI

 Manufacturer	 Serial number	 Uncoated device
 Date of manufacture	 Reference number	 Dangerous voltage
 Consult instructions for use	 Temperature limitation	 Open here
 Do not resterilize	 Lot number	 CRT-D RA, RV, LV
 Do not use if package is damaged	 Investigational use only	 ICD RA, RV
 Do not reuse	 Place telemetry wand here	 ICD RV
 Use by	 Sterilized using ethylene oxide	
  Authorized representative in the European Community		

Refer to the System/Reference Guide for all programmable parameters and values.

If the container seal is broken, Boston Scientific will not accept contents for credit or restock.

Caution: Do not store near magnets or magnetic devices.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician trained or experienced in device implant and follow-up procedures.

