

X. SIGNS AND LABELS

1. THE MOTIONPOD™

The MotionPod is the part of the system in contact with the subject.

Label of the MOTIONPOD™ (on the Black bottom cover)



15 mm of diameter

2. THE MOTIONCONTROLLER™

The label of the MOTIONCONTROLLER™ figures at the bottom of the case



Dimension~ 30 x 33 mm

XI. MAINTENANCE

The MOTIONPOD™ system does not require any particular maintenance other than cleaning after each use. Manufacturer calibration is carried out at the time of manufacture of the MOTIONPOD™. Follow instructions on section VIII for re-calibration.

The MOTIONPOD™ should not be opened. Battery replacement operation should only be carried out by supplier's qualified personnel.

FCC ID Label Location

Label here

