

ID Label / Label Location Phonak Audéo M-R

Pursuant to paragraph §2.925 of 47 C.F.R. the hearing aid family Phonak Audéo M-R, which is subject of this filing, shall bear a permanently affixed, readily visible label listing the information as specified in §2.926 and §15.19(a) of 47 C.F.R. Due to the specific use of these accessories for hearing instruments (as shown on Figure 1) their external surface is in a permanent contact with the skin behind the ear. Therefore, it is not appropriate to permanently affix a label on the outer surface of the device in order to prevent any possible skin irritations. The only possible place where the FCC ID and the serial number of the device can be put is the battery compartment. However, in this case the size of the FCC ID will be smaller than 2 points, which does not satisfy the requirements for “readily visible” label as defined in §2.925(d)(2) and §2.925(g).

Therefore, due to the very small size of the devices and pursuant to paragraph §15.19(a)(5) of 47 C.F.R., the FCC ID and the statement specified in §15.19(a)(3) are placed in the user manuals as shown in the example on Figure 3.




Figure 1: Phonak hearing aids positioning

External Pictures

Brand & Product Model labelling



Phonak Audéo M-R




Figure 2a: Sample of Phonak Audéo M-R hearing aid with a visible brand name and model number with the serial number inside of the device affixed on the special plate

Housing Colors

Phonak Audéo M-R models, which are subject of this certification, are available in the following colors



Figure 2b: colors Phonak Audéo M-R models

Figure 3: IC Certification Number, product models and compliance statement included in the user guide

Phonak Audéo M-R : English version

This user guide is valid for:

Wireless hearing aids



Phonak Audéo M90-R	2018
Phonak Audéo M90-RT	2018
Phonak Audéo M70-R	2018
Phonak Audéo M70-RT	2018
Phonak Audéo M50-R	2018
Phonak Audéo M50-RT	2018
Phonak Audéo M30-R	2018
Phonak Audéo M30-RT	2018
Phonak Audéo M-R Trial	2018
Phonak Audéo M-RT Trial	2018

Non wireless charging accessories



Phonak Charger Case Combi	2018
Phonak Power Pack	2016
Phonak Mini Charger Case	2018



The wireless models listed on page 2 are certified under:

**Phonak Audéo M-R & M-RT (M90/M70/M50/M30)
Phonak Audéo M-R Trial & M-RT Trial**

USA	FCC ID: KWC-LDR
Canada	IC: 2262A-LDR

Notice 1:

This device complies with Part 15 of the FCC Rules and with RSS-210 of Industry Canada. Operation is subject to the following two conditions:

- 1) this device may not cause harmful interference, and
- 2) this device must accept any interference received, including interference that may cause undesired operation.

Notice 2:

Changes or modifications made to this device not expressly approved by Sonova AG may void the FCC authorization to operate this device.

Notice 3:

This device has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules and ICES-003 of Industry Canada. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This device generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is

no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and receiver.
- Connect the device into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.