RF Technologies Inc.

Quick Response Nurse Call with Check In Models 0800-0539, 0800-0540, 0800-0541, and 0800-0542

FCC ID: KXU-PCFSZ242 IC: 2719A-PCFSZ242


These units are professionally installed against the wall. The installer will mount the back plate to the wall before attaching the main unit to the back plate with a screw at the bottom of the main enclosure. For long-term visibility of the compliance label — it is installed on the inside of the back plate, rather than the back side, which would be against the wall. It is a permanently visible label whenever the batteries are changed.

See Internal/External Photos for exact label placements.

Description: FCC Label, QR Premiere Pull Cord

Revision: B


Mark as follows:



Description: FCC Label, QR Premiere Pull Cord w/ Check In

Revision: B

Mark as follows:




Description: FCC Label, QR Premiere Wall Mount (red

button)

Revision: B

Mark as follows:



Description: FCC Label, QR Premiere Nurse Call Jack

Revision: B

Mark as follows:

