

## Pendant Transceiver

### Resident Confidence

The Pendant Transceiver can be used with all Quick Response® Wireless Call Solutions. The transceiver is worn around the resident's neck or clipped to a belt. This wireless device allows residents to move about your facility with the confidence that if they are in need of assistance, they can easily summon help with the press of a button.

### Continuous Assurance of Proper Operation

The Quick Response Wireless Call Solution fully supervises the transceiver and will alert staff if a problem develops or if the transceiver is outside the area of coverage. The transceiver is battery operated, and the battery can be replaced by the facility. When the battery is low, staff is continuously notified until the battery is changed.

### Advanced Direct Sequence Spread Spectrum Technology

The transceiver operates in the 2.4 GHz frequency band using Advanced Direct Sequence Spread Spectrum Technology. This technology assures reliable signal reception over a large area.



### Features

- Completely wireless
- Advanced 2.4 GHz Direct Sequence Spread Spectrum Technology
- Fully supervised
- Easy to use

### Specifications

**Power**

3.0V Lithium battery

**Battery Life**

1-2 years

**Frequency**

2.4 GHz Direct Sequence Spread Spectrum

**Frequency Range**

2.405 - 2.480 GHz

**Bandwidth**

5 MHz

**Transmit Power**

0 dBm

**Certification**

FCC, Part 15  
Zigbee 802.15.4 compliant

**Size**

2.13" H x 1.34" W x .60" D  
(5.41 cm x 3.40 cm x 12.7 cm)

**Weight**

1.4 oz. (39.69g)