

Notice - This device contains a magnetically activated switch and has been turned off before shipment to preserve the battery life. Stray magnetic fields may accidentally turn the device on. Verify device is off prior to storing (see user manual appropriate for device).

Content Notice: This device may contain some recycled electronic components. Recycled components must pass the same final inspection and testing as new components. The standard warranty is in effect on this device. PhysioTel, Data Sciences, and the Data Sciences logo are trademarks of Data Sciences International.

WARNING: FOR ANIMAL USE ONLY - NOT FOR HUMAN USE
This product may be covered by one or more of the following US patents; (7344505 and D572364)

Manufactured in the USA by Data Sciences International

P/N 006611-001 Rev05

CAUTION



DSITM

119 14th Street NW,
Suite 100, St. Paul,
MN 55112 USA
(651) 481-7400

SoHo-X00N

Implant

S/N: 1284572

P/N: 270-0700-001

Mfg. Date: 05-Jun-2024

OPID: AJH

Exp. Date: 05-Jun-2025

FCCID: MHASOHO

IC: 5681A-SOHO

CE UK CA [R] 203-JNXXXX

STERILE

S/N:

Mfg. Date: 05-Jun-2024

Model: SoHo-X00N

Implant Sensing Channels:

1 Temperature

1 Activity

Exp. Date: 05-Jun-2025

Implant Sensing Channels:

1 Temperature

1 Activity



Scan to access DSI Support

Peel Here For Your Records