

| REV | ECO | DESCRIPTION | DRFT | DATE | CHKD | APPVD |
|-----|---------|--|------|-------|----------|-------|
| 000 | ER23791 | ENG RELEASE | RB | 8/00 | | |
| A | MR26113 | MFG RELEASE | RB | 10/00 | BS | JI |
| B | 065678 | ADD CHANNEL NOS. 518A TO 520A & 600 TO 620 | RB | 10/00 | BS | MS |
| C | 065902 | ADD 'A' TO CHANL'S 600 TO 619; REM CHANL 620 | RB | 11/00 | BS | MS |
| D | 066872 | ADDED 340M SERIES STYLE | MM | 4/01 | <i>W</i> | |

CHANNEL
XXXXA

USED FOR 0340 PLUS SERIES

CHANNEL
OBXXXX


USED FOR 340M SERIES

NOTES:

- 1-PRINT COLOR BLACK.
- 2-REPLACE X'S WITH THE APPROPRIATE CHANNEL NUMBER.

PRELIMINARY

DRAWING FILENAME: 2004709-D.DWG

| | | | |
|--|-----------|--|--|
| PROPRIETARY & CONFIDENTIAL | | THIS DOCUMENT AND THE INFORMATION HEREIN ARE THE PROPERTY OF GEMS <i>IT</i> AND NEITHER MAY BE USED FOR ANY PURPOSE, NOR DISCLOSED TO ANY PARTY, WITHOUT THE WRITTEN CONSENT OF GEMS <i>IT</i> . | |
| DO NOT SCALE - WORK TO DIMENSIONS SHOWN. REMOVE BURRS, BREAK SHARP EDGES & CORNERS. | | DRAWN R. BARNES DATE 7/16/00 CHECKED B. SHOTWELL DATE 10/5/00 ENGINEER JOE IZZO DATE 10/5/00 APPROVED JOE IZZO DATE 10/5/00 | GE Medical Systems <i>Information Technologies</i>  |
| UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES. | | TITLE LABEL, CHANNEL NUMBERS TRANSMITTER & RECEIVER-340/340M TELEMETRY | |
| FRACTIONS ± N/A 2 PLACE DECIMALS ± N/A 3 PLACE DECIMALS ± N/A ANGLES ± N/A | | MATERIAL MAKE FROM 440024 | FINISH N/A |
| 2004701-001 | TELEMETRY | SIZE B | DRAWING NUMBER 2004709 |
| NEXT ASSY | USED ON | REV D | |
| APPLICATION | | SCALE 2/1 | SHEET 1 OF 1 |