

Dexcom		Filename: AW-1000422_REV001_CTN-AW_Stelo_BioSX_1PK_US_en	Agency: Cinco	Job No.: N/A
Part No./Rev.: AW-1000422_REV001		Rev. Date: 3/7/24		
Spec No./Rev.: MT-1000422_REV001		Dexcom Contact: Tiffany Carmouche		
Language Sequence: English		Department Owner: Product Labeling		
No. of Covers:	N/A <input checked="" type="checkbox"/>	Inside No. of pages:	N/A <input checked="" type="checkbox"/>	No. of sides: 1 <input checked="" type="checkbox"/> 2 <input type="checkbox"/>
Spot Colors		Process Colors		Legal/Labeling Text
Custom Spot <input checked="" type="checkbox"/> White and Plum N/A <input type="checkbox"/>		CMYK <input checked="" type="checkbox"/> Black <input type="checkbox"/> N/A <input type="checkbox"/>		Font: Arial Size: 6 pt Color: CMYK 0/0/0/70
All information contained herein is the CONFIDENTIAL property of Dexcom, Inc. and may not be duplicated or released without the expressed permission of Dexcom, Inc.				Dieline <input checked="" type="checkbox"/> N/A <input type="checkbox"/> 



**Indications for use:** The Stelo Glucose Biosensor System is an over-the-counter (OTC) integrated Continuous Glucose Monitor (CGM) intended to continuously measure, record, analyze, and display glucose levels. The Stelo Glucose Biosensor System is not a blood glucose monitoring device. The Stelo Glucose Biosensor System helps to detect normal (euglycemic) and low or high (glycemic) glucose levels. The Stelo Glucose Biosensor System may also help the user better understand how lifestyle and behavior modification, including diet and exercise, impact glucose excursion. The user is not intended to take medical action based on the device output without consultation with a qualified healthcare professional.

**Contraindications:** No MRI/CT/diathermy - MR unsafe.

**Warning:** Choking hazard. Don't use if you have problematic hypoglycemia. Don't use if you are on dialysis. Don't ignore low/high symptoms. When your sensor readings are not consistent with your symptoms, consult your healthcare provider. Consult your healthcare provider before making any medication adjustments based upon your sensor readings.

- Track glucose 24/7
- Waterproof up to 8 feet\*
- App works on Apple® and Android™ smart devices\*

Reliable glucose monitoring for long-term wellness

Turn healthy goals into healthy habits

Reveal patterns that help you turn insight into action

Dexcom, Inc.  
6340 Sequoia Drive  
San Diego, CA 92121 USA

STERILE EO

IP58

MD

FCC ID PH29788  
MT-26023-40

36 °F 100% 96 °F 10%

© 2024 Dexcom, Inc.

Glucose Biosensor Kit Contents Include

- 1 Biosensor
- 1 Applicator
- 1 Set of Instructions

UDI

Glucose Biosensor Kit (1)

by Dexcom

NOT FOR RESALE

REF

LOT

(243) MT-1000422Rev001



Reliable glucose monitoring  
for long-term wellness



Turn healthy goals into  
healthy habits



Reveal patterns that help  
you turn insight into action



→ Track glucose 24/7

→ Waterproof up to 8 feet†


→ App works on Apple® and Android™ smart devices†

**Indications for use:** The Stelo Glucose Biosensor System is an over-the-counter (OTC) integrated Continuous Glucose Monitor (iCGM) intended to continuously measure, record, analyze, and display glucose values in people 18 years and older not on insulin. The Stelo Glucose Biosensor System helps to detect normal (euglycemic) and low or high (dysglycemic) glucose levels. The Stelo Glucose Biosensor System may also help the user better understand how lifestyle and behavior modification, including diet and exercise, impact glucose excursion. The user is not intended to take medical action based on the device output without consultation with a qualified healthcare professional.

**Contraindications:** No MRI/CT/diathermy - MR unsafe.

**Warning:** Choking hazard. Don't use if you have problematic hypoglycemia. Don't use if you are on dialysis. Don't ignore low/high symptoms: When your sensor readings are not consistent with your symptoms, consult your healthcare provider. Consult your healthcare provider before making any medication adjustments based upon your sensor readings.

# stelo



Glucose Biosensor Kit (1)

by Dexcom



Dexcom, Inc.  
6340 Sequence Drive  
San Diego, CA 92121 USA



FCC ID PH29788  
MT-26023-40

If you need to track your glucose levels when your Stelo readings don't match how you feel, a fingershick may be an option as needed. †The Stelo sensor is waterproof and may be submerged under eight feet of water for up to 24 hours without failure when properly installed. ‡Smart device required and sold separately. For a list of compatible devices, visit [stelo.com/compatibility](http://stelo.com/compatibility).

Symbols glossary: [stelo.com/symbols](http://stelo.com/symbols)

Covered by patents: [stelo.com/patents](http://stelo.com/patents)

By using the Stelo sensor, you agree to the Terms of Use and Policies located at [stelo.com](http://stelo.com). To the extent allowed by applicable law, the Stelo sensor is provided to you "as is" and without any warranty, express or implied. Dexcom, Stelo, and any related logos and design marks are either registered trademarks or trademarks of Dexcom, Inc. in the United States and/or other countries.

AW-1000422 REV 001 MT-1000422

© 2024 Dexcom, Inc.

## Glucose Biosensor Kit Contents Include

- 1 Applicator
- 1 Biosensor
- 1 Overpatch
- 1 Set of Instructions

UDI

NOT FOR RESALE

REF

LOT

