

BIOTRONIK

Silicone

BioMonitor 2-AF

2015-09-30 -14°F \pm +113 °F
Use By
Stor. Temp.

12345678
Serial No.

398493
Cat. No.



(01)040354791355221171150930(2)112345678

Implantable Cardiac Monitor

ProMRI®

BioMonitor 2-AF

89
Prod ID

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner)

402-405 MHz
TP2



FCC ID: Q1BIM2

30009

Distributor:

BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
(800) 547-0394 (24-hour)
(503) 635-9936 (Fax)

Manufacturer:

BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany
Tel +49 (0) 30 68905-0
Fax +49 (0) 30 6852804
sales@biotronik.com

BIOTRONIK

The device and the included tool(s) are intended for single use only, have been sterilized with ethylene oxide and are not resterilizable. Refer to Technical Manual for instructions.



BioMonitor 2-AF

398493
Cat. No.

12345678
Serial No.

2015-09-30
Use By

89
Prod ID



BIOTRONIK

BioMonitor 2-AF
SN 12345678
PID 89
BIOTRONIK

BioMonitor 2-AF
SN 12345678
PID 89
BIOTRONIK

BioMonitor 2-AF
SN 12345678
PID 89
BIOTRONIK

BioMonitor 2-AF
SN 12345678
PID 89
BIOTRONIK

BioMonitor 2-AF
SN 12345678
PID 89
BIOTRONIK

BioMonitor 2-AF
SN 12345678
PID 89
BIOTRONIK

BioMonitor 2-AF



Implantable Cardiac Monitor

ProMRI®

398493
Cat. No.


2015-09-30
Use By

12345678
Serial No.

89
Prod ID

STERILE

Sterilized with ethylene oxide



For single use only, do not resterilize. Do not use if packaging is damaged. Refer to technical manual for instructions.

BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany

BIOTRONIK



(01)040354791355221171150930(2)112345678

12345678
SN
89
PID

12345678
SN
89
PID



AD13FF

BioMonitor 2-AF

BIOTRONIK



Implantable Cardiac Monitor

398493
Cat. No.

12345678
Serial No.

2015-09-30
Use By

89
Prod ID

ProMRI®

Individual Data Sheet



Parameter Values at Factory Release

Parameter	Values at Factory Release
Atrial fibrillation detection	OFF
High ventricular rate detection	OFF
Bradycardia detection	OFF
Sudden rate drop detection	OFF
Asystole detection	OFF
Patient trigger	OFF
Home monitoring	OFF

For single use only. Refer to Technical Manual for instructions.

Distributor:
BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
(800) 547-0394 (24-hour)
(503) 635-9936 (Fax)

BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany
Tel +49 (0) 30 68905-0
Fax +49 (0) 30 6852804
sales@biotronik.com

BIOTRONIK

BIOTRONIK