


# Lumax 640 DR-T



IS-1  
IS-1  
DF-1  
DF-1

REF 381470  
Cat. No.

SN 12345678  
Serial No.

PID 89  
Prod.ID

Use By 2012-12

Store Temp. 41-113 °F





(01)04035479122874(17)121231(21)12345678

Implantable Cardioverter Defibrillator

# Lumax 640 DR-T

ProMRI®


Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licensed practitioner)



Manufacturer:  
BIOTRONIK SE & Co. KG  
Weermannkehrle 1  
12359 Berlin - Germany  
Tel +49 (0) 30 68905-0  
Fax +49 (0) 30 6852804  
sales@biotronik.com


Distributor:  
BIOTRONIK, Inc.  
6024 Jean Road  
Lake Oswego, OR 97035-5369  
(800) 547-0394 (24-hour)  
(503) 635-9936 (Fax)

FCC ID: ORLUMAXT50



For single use only, do not resterilize. Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

# Lumax 640 DR-T



IS-1  
IS-1  
DF-1  
DF-1

Implantable Cardioverter Defibrillator

ProMRI®


REF 381470  
Cat. No.

Use By 2012-12

SN 12345678  
Serial No.

PID 89  
Prod.ID

STERILE EO  
Sterilized with ethylene oxide



For single use only, do not resterilize. Do not use if packaging is damaged. Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

BIOTRONIK SE & Co. KG  
Weermannkehrle 1  
12359 Berlin - Germany




(01)04035479122874(17)121231(21)12345678

SN 12345678

SN 12345678

# Lumax 640 DR-T



IS-1  
IS-1  
DF-1  
DF-1

Implantable Cardioverter Defibrillator

Use By 2012-12  
Prod.ID 89

REF 381470  
Cat. No.  
SN 12345678  
Serial No.

# Individual Data Sheet



Parameter	Values at Factory Release	Programmable Parameter Ranges
Bradycardia Therapy Mode	OFF	AAI(R); V00; VVI(R); VDI(R); VDD(R); D00; DDI(R); DDD(R); OFF
Basic Rate; min-1		30...(5)...100...(10)...160
Pulse Amplitude; V		0.5...(0.25)...4...(0.5)...6; 7.5
Pulse Width; ms		0.4; 0.5...(0.25)...1.5
ATPs	OFF	OFF; BURST; RAMP; SCAN
Shock	OFF	OFF; DEFIBRILLATION with/without Confirmation
X-Ray Identification	RH	

For single use only. Refer to Technical Manual for Instructions.

Distributor:  
BIOTRONIK, Inc.  
6024 Jean Road  
Lake Oswego, OR 97035-5369  
(800) 547-0394 (24-hour)  
(503) 635-9936 (Fax)

BIOTRONIK SE & Co. KG  
Weermannkehrle 1  
12359 Berlin - Germany  
Tel +49 (0) 30 68905-0  
Fax +49 (0) 30 6852804  
sales@biotronik.com