

41-113 °F  
 Ster. Temp.

2013-12  
 Use By

SN 12345678  
 Serial No.

REF 391917  
 Cat. No.




(01)104035479 128586(17)1131231(21)112345678

ProMRI®



**Iforia 3 DR-T**

PID 89  
 Prod.ID




Implantable Cardioverter Defibrillator

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician for properly licensed practitioner)



FCC ID: QRTACHNXT



For single use only, do not resterilize. Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

Manufacturer:  
 BIOTRONIK SE & Co. KG  
 Woermannkehe 1  
 12359 Berlin - Germany  
 Tel. +49 (0) 30 68905-0  
 Fax +49 (0) 30 6852804  
 sales@biotronik.com

Distributor:  
 BIOTRONIK, Inc.  
 6024 Jean Road  
 Lake Oswego, OR 97035-5369  
 (800) 547-0394 (24-hour)  
 (503) 635-9936 (Fax)

**Iforia 3 DR-T**

REF 391917

Cat. No.

SN 12345678

Serial No.

2013-12

Use By

PID 89

Prod.ID



Iforia 3 DR-T  
 SN 12345678  
 PID 89  
 BIOTRONIK

Iforia 3 DR-T  
 SN 12345678  
 PID 89  
 BIOTRONIK

Iforia 3 DR-T  
 SN 12345678  
 PID 89  
 BIOTRONIK

Iforia 3 DR-T  
 SN 12345678  
 PID 89  
 BIOTRONIK

Iforia 3 DR-T  
 SN 12345678  
 PID 89  
 BIOTRONIK

Iforia 3 DR-T  
 SN 12345678  
 PID 89  
 BIOTRONIK

Iforia 3 DR-T  
 SN 12345678  
 PID 89  
 BIOTRONIK

**Iforia 3 DR-T**

IS-1  
 IS-1  
 DF-1  
 DF-1



ProMRI®

Implantable Cardioverter Defibrillator

REF 391917

Cat. No.

SN 12345678

Serial No.

2013-12

Use By

PID 89

Prod.ID

STERILE EO

Sterilized with ethylene oxide



For single use only, do not resterilize. Do not use if packaging is damaged. Refer to technical manual for instructions. ICD is shipped with all therapy disabled.

BIOTRONIK SE & Co. KG  
 Woermannkehe 1  
 12359 Berlin - Germany



(01)104035479 128586(17)1131231(21)112345678

SN 12345678

SN 12345678

**Iforia 3 DR-T**

IS-1  
 IS-1  
 DF-1  
 DF-1



Implantable Cardioverter Defibrillator

2013-12

Use By

PID 89

Prod.ID

REF 391917

Cat. No.

SN 12345678

Serial No.



**Individual Data Sheet**

Parameter	Values at Factory Release	Programmable Parameter Ranges
Bradycardia Therapy Mode	OFF	AAIRI; V00; VVIRI; VDIRI; VDDR; D00; DDIRI; DDDIRI; OFF
Basic Rate; min-1		30...151...100...101...160
Pulse Amplitude; V		0.5...10.251...4...10.51...6; 7.5
Pulse Width; ms		0.4; 0.5...10.251...1.5
ATPs	OFF	OFF; BURST; RAMP; SCAN
Shock	OFF	OFF; DEFIBRILLATION with/without Confirmation
X-Ray Identification	NT	

For single use only. Refer to Technical Manual for Instructions.

Distributor:

BIOTRONIK, Inc.  
 6024 Jean Road  
 Lake Oswego, OR 97035-5369  
 (800) 547-0394 (24-hour)  
 (503) 635-9936 (Fax)

BIOTRONIK SE & Co. KG  
 Woermannkehe 1  
 12359 Berlin - Germany  
 Tel. +49 (0) 30 68905-0  
 Fax +49 (0) 30 6852804  
 sales@biotronik.com

