

Shipping box Label - RIASJMRFB



Cardiac Rhythm Management Division | 15900 Valley View Court | Sylmar, CA 91342 USA | +1 818 362 6822 | 800 722 3774 (US Only)

SN
Serial number:



Use before:

ACCENT™ RF DR
Pulse Generator

PM2210

DDDR
IS-1



ACCENT™ RF DR Pulse Generator

PM2210



Use before:

SHIPPED SETTINGS

Mode	Rate	AV Delay	Max Track	Sensor
DDD	60 ppm	200ms	110 ppm	Off

	ATRIAL	VENTRICULAR
Pace/Sense	Bipolar/Bipolar	Bipolar/Bipolar
Pulse Amplitude	2.5 V	2.5 V
Pulse Width	0.4 ms	0.4 ms
Sensitivity	0.5 mV	2.0 mV
Refractory	275 ms	250 ms

CONTENTS

One pulse generator, accessories and product literature.

STERILE EO

Contents of this package have been ethylene oxide sterilized. Sterility cannot be guaranteed if sterile package is damaged or opened. Refer to the accompanying manual.

STORAGE CONDITIONS

Store at -5°C (23°F) to 50°C (122°F).




CAUTION

Rx only

Proposition 65, a State of California voter initiative, requires the following notice:

WARNING: This product and its packaging have been sterilized with ethylene oxide. This packaging may expose you to ethylene oxide, a chemical known to the state of California to cause cancer or birth defects or other reproductive harm.

DDDR
IS-1



Serial number:

SN

Manufacturing date:



FCC ID: RIASJMRFB
IC: 7067A-SJMRF2

Single Use Only!

ACCENT™ RF DR
Pulse Generator

PM2210



Use before:

DDDR
IS-1

Serial number:

SN

Note: The shipping box label is configured the same for all Accent and Anthem RF implant devices.

Label location- RIASJMRFB

The FCC ID is located on the shipping box label for the RF enabled implantable devices.

