

Labeling on device

For US

The back of the UPD-3 is shown in figure 1 below. FCC Identifier will be placed where figure 1 indicate “Affix FCC Label”. The detail is shown in figure 2.

For Canada

The back of the UPD-3 is shown in figure 1 below. IC Identifier will be placed where figure 1 indicate “Affix IC Label”. The detail is shown in in figure 2.




Figure 1 Location of the label

**ENDOSCOPE POSITION
DETECTING UNIT
MODEL UPD-3
POWER INPUT 100-240V ~
50/60Hz 110VA**

Height of font

3mm

TV TYPE NTSC
FCC ID:S8QUPD-3
IC:4763B-UPD3

OLYMPUS MEDICAL SYSTEMS CORP.
2951 Ishikawa-cho, Hachioji-shi, Tokyo, Japan
MADE IN JAPAN

Figure 2 Datal of the label