

# EDAN Patient Monitor

Model: iM60

 YYYY-MM

 0123

SN



CCCC

Anti-electro-shock Type: Class I

100V-240V~ 50Hz/60Hz 1.4A-0.7A

FCC ID: SMQIM60EDAN



EC REP

Shanghai International Holding Corp.GmbH(Europe)

Eiffestrasse 80,D-20537 Hamburg Germany

Tel:+49-40-2513175 Fax:+49-40-255726

IPX1

Rx only (U.S.)



Edan Instruments, Inc.

3/F - B, Nanshan Medical Equipments Park, Nanshan

Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China

Tel: +86-755-26898326 Fax: +86-755-26898330

[http:// www.edan.com.cn](http://www.edan.com.cn)

## Label and Label Location

FCC ID : SMQIM60EDAN

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

### Label Location



## Label and Label Location

Label:

FCC ID : SMQIM60EDAN

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:  
(1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Label Location

