

EXHIBIT A - FCC ID LABEL INFORMATION

FCC ID Label & Label Location

GROUNDING: Reliability can only be achieved when connected to an equivalent receptacle marked **HOSPITAL GRADE**.


WARNING: Disconnect AC power cord before removing the battery cover.


CAUTION: Electrical shock hazard. Do not open case. Refer servicing to qualified personnel.


CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.

DANGER-Possible explosion hazard if used in the presence of flammable anesthetics.


DANGER-Risque d'explosion possible si employe en presence d'anesthésiques inflammables.

 UL 60061-1
CSA 901.1
MCN 180882

 Type CF

 EN-2

IPX1
Medical Equipment
Patents Pending

 FCC ID: 513-20708
IC: 5827A-20708

**120 VAC, 50-60 Hz,
50 VA Max**

FDA regulation 21 CFR 821 requires that this device be tracked when placed in home care use. The home care distributor must register this device with Hospira.

Hospira, Inc.,
Lake Forest, IL 60045 USA 735-54826-001

Rear View of EUT/FCC ID Label

