



GROUNDING: Reliability can only be achieved when connected to an equivalent receptacle marked **HOSPITAL GRADE**.

WARNING: Disconnect AC power cord before removing the battery cover.

CAUTION: Electrical shock hazard. Do not open case. Refer servicing to qualified personnel.

CAUTION: Federal (USA) law restricts this device to sale by or on the order of a physician or other licensed practitioner.

DANGER-Possible explosion hazard if used in the presence of flammable anesthetics.

DANGER-Risque d'explosion possible si employe en presence d'anesthésiques inflammables.

UL 60601-1
 @ CSA 601.1
 MCN 160992
 c US


 Type CF

✓ EN-2

IPX1
 Medical Equipment
 Patents Pending

This Device Contains
 FCC ID: STJ80411396001
 IC: 5627A-80411396
 Model: CUSTOM DWL-AG132



120 VAC, 50-60 Hz,
 50 VA Max



FDA regulation 21 CFR 821 requires that this device be tracked when placed in home care use. The home care distributor must register this device with Hospira.

Hospira, Inc.,
 Lake Forest, IL 60045 USA 735-04828-001, A