

Angel Medical Systems, Inc.

**Guardian® IMD/EXD System for monitoring electrograms (EGMs),
402 MHz – 405 MHz**

FCC ID: THL- IMDAG101

Photo One – External

Implantable Medical Device (IMD)

Sealed Product with Measurement / Antenna Leads

FCC ID Number: Printed on Sterile Package of the IMD Implant.

Once IMD is Implanted, FCC ID not visible.



Angel Medical Systems, Inc.

**Guardian[®] IMD/EXD System for monitoring electrograms (EGMs),
402 MHz – 405 MHz**

FCC ID: THL- IMDAG101

Photo Two – External

Implantable Medical Device (IMD)

Sealed Product / Backside with out Leads

**(Note: “Not for human use” labeling is used for verification type tests,
not on final product.**



Angel Medical Systems, Inc.

**Guardian[®] IMD/EXD System for monitoring electrograms (EGMs),
402 MHz – 405 MHz**

FCC ID: THL- IMDAG101

Photo Three – External

**Implantable Medical Device (IMD)
External Antenna Connection – Close-up**



Angel Medical Systems, Inc.

Guardian[®] IMD/EXD System for monitoring electrograms (EGMs), 402 MHz – 405 MHz

FCC ID: THL- IMDAG101

Photo Four – External

Implantable Medical Device (IMD)

Printed Art Work of IMD Implant Sterile Packaging
with FCC ID Number ... THL-IMDAG101

Note: The IMD is a Sterile Product, implanted and not visible during operation.

AngelMed Guardian[®] IMD Model # AG101



Implantable Medical Device (IMD) for detecting and warning a patient of a variety of changes in his or her electrogram.

For use ONLY with Angel Medical Systems, Inc. Model LA-001 Lead Adapter and St. Jude Model 1488T Bipolar Cardiac Lead

Contents: One IMD, One Torque Wrench

The contents (STERILE) of this inner package have been sterilized by ethylene-oxide gas. Sterility cannot be guaranteed if sterile package is damaged or opened.



STERILE EO



YYYY-MM



YYYY-MM



SN

XXXXXX



-10°C/14° F and +55°C/131°F

Caution: Investigational Device. Limited by federal (or United States) law to investigational use.

Manufactured for: Angel Medical Systems, Inc.
1 Sheila Drive
Tinton Falls, NJ 07724 USA
Phone: 732-212-1888 (USA)



THL-000AG101

(Hi-Tronics Document No.)